

# HPLC

## Column Protection Guide

Version 0113

### **Includes:**

- Mobile Phase Limitations
- Column Storage Tips
- Column Protection Devices

# COLUMN PROTECTION GUIDE

## TABLE OF CONTENTS

### **INTRODUCTION TO COLUMN CARE AND USAGE 1-4**

Selecting the Right Tubing and Fittings  
Column Installation

### **PART I - SILICA-BASED & TWIN™ TECHNOLOGY COLUMNS 4-7**

Running Parameters  
Mobile Phase Considerations  
Stationary Phase Considerations  
Backpressure and Flow Rates

Column Equilibration  
Column Storage  
Column Cleaning Procedures

### **PART II - ONYX SILICA MONOLITH COLUMNS 7-10**

Running Parameters  
Mobile Phase Considerations  
Stationary Phase Considerations  
Hardware Considerations  
Column Installation

Column Equilibration  
Cleaning & Regeneration Procedure  
Column Storage

### **PART III - SFC (SUPERCRITICAL FLUID CHROMATOGRAPHY) 11**

Running Parameters  
Equilibrating Column  
Mobile Phase Considerations

Cleaning Procedure  
Column Storage

### **PART IV - AXIA PREPARATIVE COLUMNS 11-12**

Running Parameters  
Mobile Phase Considerations

Cleaning Procedure  
Column Storage

### **PART V - LUX CHIRAL COLUMNS 12-13**

Running Parameters  
Mobile Phase Considerations

Extending Lifetime and Reconditioning  
Column Storage

### **PART VI - CHIREX CHIRAL COLUMNS 13-14**

Running Parameters  
Mobile Phase Considerations

### **PART VII- BIOSEP AND YARRA SEC COLUMNS 14**

Running Parameters  
Mobile Phase Considerations  
Cleaning Procedure  
Column Storage

### **PART VIII- REZEX POLYMER-BASED COLUMNS 14-19**

Running Parameters  
Mobile Phase Considerations  
Cleaning Procedure  
Table of Specifications and Operating Parameters

### **PART IX - POLYSEP-GFC-P COLUMNS 20**

Running Parameters  
Mobile Phase Considerations  
Cleaning Procedure  
Column Storage

### **PART X - PHENOGEL GPC COLUMNS 21, 24-25**

Specifications  
Sample Considerations  
Column Storage

Solvent Switching Considerations  
Solvent Compatibility Chart

### **SOLVENT MISCIBILITY TABLE 22-23**

### **PART XI - POLYMERX RP COLUMNS 26**

Specifications  
Running Parameters  
Mobile Phase Considerations

Cleaning Procedure  
Column Storage

### **PART XII - HPLC COLUMN PROTECTION AND PERFORMANCE TESTING 26-37**

### **PART XIII - VIALS 38**

### **PART XIV - SOLID PHASE EXTRACTION (SPE) 40**

### **PART XV - HPLC ACCESSORIES 41**

### **PHENOMENEX WARRANTY 42**

# INTRODUCTION

Every Phenomenex HPLC column is a precision product which, though delicate, will provide excellent performance, reproducibility and column lifetime if cared for properly. The information and recommendations contained in this manual are designed to guide you in the care and use of your column, but should not be considered absolute. Please follow the instructions herein to maximize column performance and lifetime. Should you have any questions, please contact your Phenomenex Technical Representative or local distributor.


## UPON RECEIPT OF THE COLUMN

- **Verify the column you received is the column you ordered**
- **Check the column for physical damage which may have occurred during shipping**
- **Test the column immediately to verify performance and quality**
- **All columns are shipped in the testing solvent, unless otherwise specified**

Each Phenomenex manufactured HPLC column is individually packed and tested to ensure high column quality. Every column is supplied with its Test Chromatogram and a Specification Sheet which indicates column serial number and identity, testing conditions and operating parameters.

The warranty period begins upon receipt of the column. Testing is especially important if the column is to be placed in storage. Test the column using the same conditions in the test chromatogram. Use the formulae in Figure 1 to determine column efficiency and peak asymmetry.

Chromatographic performance depends on the entire system, not just the column. Columns are QC tested using optimum conditions to minimize bandspreading from "extra column effects." Most variations from the Phenomenex test data are due to extra-column effects created by the design of your system (i.e., injector, flow cell, connecting tubing, etc.). If you have any questions regarding your test results or the column quality, or if there are signs of damage, CONTACT PHENOMENEX OR YOUR LOCAL DISTRIBUTOR IMMEDIATELY.



*Formulae for calculating efficiency and peak asymmetry*




## SELECTING THE RIGHT TUBING AND FITTINGS

The tubing and fittings on an HPLC system contribute to system dead volume. If not minimized, dead volume can lead to band broadening and peak degradation. Please use the following guideline to keep system dead volume to a minimum and to help ensure optimum column performance.



### TUBING

The choice of tubing material is based on its chemical resistivity, application and HPLC system considerations (i.e. flow rate, backpressure, etc). Please refer to Tables 1-3 for specifics.

### TUBING DIAMETER

TABLE 1	
<b>High Pressure Tubing:</b>  1/16 in. OD x 0.010 in. ID	
<b>Inlet/Outlet Low Pressure Tubing:</b>	
 1/16 in. OD x 0.030 in. ID	 1/8 in. OD x 0.062 in. ID

### TUBING COMPATIBILITY

TABLE 2		
<b>Stainless Steel</b> (Type 316)		AVOID high concentrations of acids or halogenated salts
<b>PEEK</b> (biocompatible)		AVOID 100 % THF, chlorinated solvents, high concentrations of acids
<b>Titanium</b> (biocompatible)		Compatible with nearly all chemicals


### TUBING APPLICATIONS

TABLE 3		
Tubing ID (Inch)	Column IDs (mm)	Typical Flow Rates (mL/min)
0.002	0.30 (Fused Silica)	0.001 - 0.02
0.005	1.0 (Stainless Steel)	0.02 - 0.1
0.007	2.0 - 4.6	0.2 - 2.0
0.010	3.2 - 7.8	0.5 - 5.0
0.020	10.0 - 21.2	2.0 - 50.0
0.040	21.2 - 100.0	10.0 - 200.0

## FITTINGS

All Phenomenex column end-fittings are female inverted (internal type) with 10-32 type threading:


- The end-fitting can fit any  $\frac{1}{16}$  in. OD tubing (see page 2 for tubing considerations)
- A 10-32 threaded male nut and ferrule or a polymeric fingertight male nut is used to swage or tighten the tubing onto the fitting (see Figure 2)



## INSTALLATION CONSIDERATIONS:

- The shape of the swaged ferrule can differ between manufacturers. For Phenomenex columns, you may use Phenomenex or Valco type ferrules.
- **VERY IMPORTANT:** The seating depth of the nib (Figure 2) for Phenomenex columns is 0.080 in. Tubing **MUST** be seated all the way down into the column end-fitting. Failure to do so will result in having a small mixing chamber at the top or bottom of the column. This will lead to degraded chromatography.

\*Polymeric fingertight fittings are easy to use. They come in one piece, require **NO** tools for attachment and easily conform to the shape of the column end-fitting.



---

## COLUMN INSTALLATION

***IT IS HIGHLY RECOMMENDED THAT YOU READ THIS GUIDE FOR SPECIFIC COLUMN CONSIDERATIONS BEFORE PROCEEDING WITH THE INSTALLATION (PARTS I-XI)***

- Flush HPLC pump and line thoroughly with filtered and degassed mobile phase (without any buffers). Make sure there are no air bubbles in the system.
- Connect the column to the injector corresponding to the direction of the flow label (located on the column). Leave the outlet of the column unattached.
- Set pump to flow at 0.1 mL/min (or lowest setting) and increase to normal flow rate over 5 minutes.
- Stop flow when there is a free flow of solvent from the column outlet, wipe the end and attach to the detector
- Equilibrate the column by passing approximately 10-30 column volumes of mobile phase at normal flow rate.
- **For those columns that can be used under reversed-phase or normal phase conditions (i.e., -CN or -NH<sub>2</sub>), flush with 20-30 column volumes of IPA or THF as the intermediate solvent when switching from reversed-phase to normal phase modes, or vice versa.**

## PART I - SILICA-BASED & TWIN™ TECHNOLOGY COLUMNS

### RUNNING PARAMETERS

- Unless otherwise specified, all porous silica and TWIN technology columns should be limited to backpressure below 3500 psi (245 bar). For applications where such pressures may be exceeded, please consult Phenomenex technical support for method specific considerations. Guidelines presented within this section also apply to High Speed Technology (HST) silica-based columns. See the core-shell column guide for specific UHPLC operational tips.
- Avoid any sudden pressure changes
- If high backpressure is observed, reverse flush the column (do not try this on other manufacturers' columns)
- Use a backpressure regulator if you are experiencing out-gassing problems in the detector cell.
- Unless otherwise specified, maximum operating temperature for all silica columns is 60 °C. For core-shell and polymer specific columns, see specific guides or contact Phenomenex technical support for method specific considerations.

---

### MOBILE PHASE CONSIDERATIONS

- Use only HPLC grade solvents and water
- Use only highest purity chemicals and reagents
- Filter and degas all mobile phases prior to use
- **Make sure solvents are miscible**

Trace impurities can dramatically degrade HPLC columns. When changing to a different mobile phase, make sure the solvents and/or buffers are miscible (see Table 11). Using solvents that are immiscible with the solvent in the column can permanently damage the column. **Salt and buffer precipitation from the mobile phase can permanently damage the column.** Always check sample solubility and if possible use the mobile phase as the diluent (sample solvent).

---

### STATIONARY PHASE CONSIDERATIONS

- Maintain pH between 2.0 and 8.0\*
- Use presaturator columns and guard columns
- Avoid aldehydes and ketones with amino columns

Silica-based columns are pH sensitive. Low pH ( $\leq 2.0$ ) will hydrolyze the bonded phase (strip off the functional groups) and high pH ( $\geq 8.0$ ) will dissolve the silica. If the mobile phase pH is near 2.0 or 8.0, use a presaturator column.

\*Consult Phenomenex for columns that have extended pH ranges.

## BACKPRESSURE AND FLOW RATES

The following backpressures are typically observed during QC testing of the following particle sizes and dimensions.

TABLE 4

Particle Size $\mu\text{m}$	Internal Diameter(mm)	Typical Flow Rate(mL/min)	Typical Pressure(psi)	
			150 mm*	250 mm*
1.7	2.1	0.3	6700	NA
2.6	2.1	0.2	6400	NA
2.6	3.0	0.8	5500	NA
2.6	4.6	1.85	5000	NA
3	2.0	0.2	1500	2400
3	3.0	0.3	1500	2400
3	4.6	0.75	1500	2300
5	2.0	0.2	650	1000
5	3.0	0.5	900	1400
5	4.6	1.0	850	1200
10	10.0	5.0	900	1000
Axia Luna	21.2	20.0	350	500

\* column length

Columns can be operated at any flow rate that is consistent with the backpressure limitations described previously. Flow rates should be optimized to provide the best efficiency for your sample.

## SCALING UP/SCALING DOWN

Adjusting flow rates for different column internal diameters is straightforward. To keep the retention times constant, the flow rates and loading capacity must be adjusted according to the column's internal diameter. Assuming column length does not change:

$$X = \text{Scale Factor} = \frac{(\text{radius column B})^2}{(\text{radius column A})^2}$$

From a 4.6 mm ID column some approximate scaling factors are:

TABLE 5

Internal Diameter	Scaling Factor
1.0 mm	0.05x
2.0 mm	0.2x
3.0 mm	0.5x
10.0 mm	5x
21.2 mm	21x



HPLC columns running water-free, flammable organic solvents (e.g., normal phase, chiral, GPC) can generate static electricity and should be properly grounded to avoid a potentially dangerous electrical discharge.

## COLUMN STORAGE

- Column storage conditions affect column lifetime
- Never store columns containing buffers or ion-pairing reagents
- Flush with 5 column volumes of mobile phase without buffer to remove any buffers or salts

Storage Conditions for Silica-Based HPLC Columns:

TABLE 6

Column Type	Storage Solvent
Reversed Phase C18, C12, C8, C4, C2, C1, Phenyl, PFP	65 % Acetonitrile/ 35 % Water
Normal Phase Silica, CN, NH <sub>2</sub> , PAC, Diol Alumina	Isopropanol or Hexane
Ion-Exchange SAX, SCX, WAX, WCX	Methanol*
Size-Exclusion Diol	0.05 % NaN <sub>3</sub> in water or 10 % methanol
HILIC Luna HILIC	80 % Acetonitrile/ 20 % Water

\*Flush column with 50 mL HPLC grade water prior to storage solvent

---

## COLUMN CLEANING PROCEDURES

The following conditions apply to Phenomenex silica-based columns with the exception of chiral columns (see Parts V and VI):

- Before starting any kind of cleaning procedure, make sure your in-column solvent or mobile phase is miscible with the recommended cleaning solvent(s).
- Flow rates should be  $1/5 - 1/2$  of the typical flow rate.

To estimate the column volume, use the following

equation:  $V = \pi r^2 L$

$V$  = column volume in mL

$r$  = column radius in cm

$L$  = column length in cm

### UNBONDED SILICA (Si)

Rinse with 10 Column Volumes each of:

- Hexane
- Methylene Chloride
- Isopropanol
- Methylene Chloride
- Mobile Phase

### Water Removal Procedure:

Flush column with 30 mL 2.5 % 2,2-dimethoxy-propane and 2.5 % glacial acetic acid in Hexane

### BONDED NORMAL PHASE (CN, NH<sub>2</sub>, DIOL, PAC)

Rinse with 10 Column Volumes each of:

- Chloroform
- Isopropanol
- Methylene Chloride
- Mobile Phase

*Exception: Luna Amino in reversed phase mode.*

### HILIC

Rinse with 10 Column Volumes each of:

- 95 % Water/5 % Acetonitrile (for buffer removal)
- 95 % 100 mM Ammonium Acetate, pH 5.8/5 % Acetonitrile
- 95 % Water/5 % Acetonitrile
- Mobile Phase

### REVERSED PHASE

**(C18, C12, C8, C5, C4, C2, C1, PHENYL, PFP, CN, NH<sub>2</sub>)**

Rinse with 10 Column Volumes each of:

- 95 % Water/5 % Acetonitrile (for buffer removal)
- THF
- 95 % Acetonitrile/5 % Water
- Mobile Phase

### REVERSED PHASE PROTEIN/PEPTIDE

**(C18, C12, C8, C5, C4, Phenyl)**

Rinse with 20 Column Volumes of mobile phase with buffer removed. Run gradient (2x):

- A) 0.1 % TFA in water
- B) 0.1 % TFA in Acetonitrile/Isopropanol (1:2)  
25 % B to 100 % B for 30 minutes

Equilibrate with 10 column volumes of mobile phase

Do not store column in TFA

### ION-EXCHANGE (SAX, SCX, WAX, WCX)

Rinse with 10 Column Volumes each of:

- 500 mM Phosphate Buffer pH 7
- 10 % Acetic Acid (Aq)
- 5 Column Volumes of Water
- 10 Column Volumes of Phosphate Buffer pH 7
- 5 Column Volumes of Water
- 10 Column Volumes of Methanol
- 10 Column Volumes of Water

### For Protein Removal

Follow the above procedure with this exception:

Substitute 10 Column Volumes of Methanol with 10 Column Volumes of 5 M Urea **or** 5 M Guanidine Thiocyanate



## NH<sub>2</sub> for HILIC or ION-EXCHANGE

Rinse with at least 20 Column Volumes each of:

- 50/50 Organic (Acetonitrile or Methanol) / 20 mM ammonium bicarbonate pH 10 (to clean ionically-bound components)
- Water
- Mobile Phase

This column cleaning procedure should only be done infrequently, as repeated exposure to high pH solutions will cause silica dissolution resulting in peak shape issues.

## GFC/SEC (Yarra SEC, BioSep SEC\*) *\*See Part VII for more details*

Rinse with 5 column volumes of 0.1 M Phosphate Buffer pH 3.0. For strongly retained proteins, run the following gradient: 100 % Water to 100 % Acetonitrile to 100 % Water over 60 minutes OR wash with 5 column volumes of 6 M Guanidine Thiocyanate or 10 % DMSO. **Do not backflush column.**

---

## SAMPLE CONSIDERATIONS

Always prefilter samples with Phenex™ syringe filters to avoid particulate contaminants which may clog the column. Use of a guard column is highly recommended to prolong the life of your analytical to preparative column.

# PART II - ONYX SILICA MONOLITH COLUMNS

## RUNNING PARAMETERS

### Pressure

- Keep backpressures below 3000 psi (200 bar)
- Avoid any sudden pressure changes
- Use a backpressure regulator if you are experiencing out-gassing problems in the detector cell
- If high backpressure is observed, reverse flush the column

Reverse the flow periodically to prevent particles and non-eluting sample components from accumulating on the column. When reversing the flow, flush the column before connecting it to the detector.

### Temperature

- Maximum operating temperature is 45 °C

As with particulate columns, preheating the mobile phase to the same temperature as the column is recommended. This can be done by placing the connecting tubing inside the column oven. [R.G. Wolcotte et al. (2000), J. Chromatogr. A., 869, 211-230].

## SPECIFICALLY for 150 x 0.1 mm dimension

### Pressure

- < 300 bar

### Flow Rate

- We recommend 1–3 µL per minute to maximize column performance

## SPECIFICALLY for SemiPrep (10 mm ID) dimension

### Pressure

- Maximum operating pressure is 150 bar (2175 psi)
- When switching valves are used – Maximum operating pressure is 100 bar (1450 psi) due to pressure fluctuations (or pressure spikes) that occur when using these

### Flow Rate

- 5 – 35 mL/min
- Fast flow rates require fast system settings.

If sampling valve or fraction collector valves switch slowly, the high flow rate of mobile phase will cause large pressure build-up momentarily. This may damage the column.

---

## MOBILE PHASE CONSIDERATIONS

- Use only HPLC grade solvents and water
- Use only highest purity chemicals and reagents
- Degas and filter mobile phases prior to use
- Make sure solvents are miscible
- For best performance use acetonitrile/water mobile phases

### Organic Solvents

Onyx columns can be used with all commonly used HPLC grade organic solvents with the following restrictions. The mobile phase should NOT contain more than 50 % tetrahydrofuran (THF), 5 % Dichloromethane (DCM), or 5 % Dimethylsulfoxide (DMSO).

Pure DMSO (up to 100  $\mu$ L) can be used as a solvent for samples. For DMSO injection volumes larger than 100  $\mu$ L, we recommend using a mixture of 50 % DMSO and 50 % diluting solvent (i.e. – Methanol).

### Buffers, Organic Modifiers, & Ion-Pair Reagents

Buffers, organic modifiers, and ion-pair reagents present no problems as long as the appropriate pH range is not exceeded. Ion-pair reagents are often difficult to completely flush from the column. Therefore, columns used with these reagents should be dedicated to the particular analysis involved.

### Acids & Bases

Do not use strong acids (i.e., hydrochloric, nitric, and sulfuric acids) in the column. Limit your use of strong bases (i.e., sodium, potassium, ammonium hydroxide) to amounts needed to adjust the pH of the mobile phase.

---

## STATIONARY PHASE CONSIDERATIONS

- Maintain pH between 2.0 and 7.5\*
- Use guard columns

\*Silica-based columns are pH sensitive. Low pH ( $\leq 2.0$ ) will hydrolyze the bonded phase (strip off functional groups) and high pH ( $\geq 8.0$ ) will dissolve the silica. If the mobile phase pH is near 2.0 or 8.0, use a presaturator column.

---

## HARDWARE CONSIDERATIONS

Onyx columns are clad with a PEEK polymer. The endfittings are also made of PEEK. DO NOT remove the endfittings from the column.

---

## COLUMN INSTALLATION

### SPECIFICALLY for 150 x 0.1 mm dimension

The tubing and fittings contribute to system dead volume. Therefore dead volume should be minimized. Onyx capillary column packages are equipped with PEEK  $1/16$  inch fittings and green sleeves. The fittings and sleeves fit with any 360  $\mu$ m OD fused silica tubing.

### Connection to Injector

- Always use in the flow direction indicated by the arrow on the column label
- Connect a 360  $\mu$ m OD fused silica tubing with the PEEK  $1/16$  inch fittings and the green sleeves to tighten the tubing

- Make sure that the tubing is seated all the way down into the fittings
- Make sure that the tubing enters all the way to the bottom of the injector port
- Keep the tubing as short as possible to avoid dead volumes
- Be careful not to over bend the capillary

### **Connection to the Detector**

Onyx capillary columns can be directly connected to any nano/capillary-HPLC UV detector (e.g. equipped with a nano/capillary flow cell) or MS. Connections are made with either a fingertight fitting or with PTFE tubing.

## **COLUMN EQUILIBRATION**

Columns can dry out during shipping and stocking, therefore thoroughly activate the Onyx packing material by equilibrating your column.

Onyx columns are shipped in acetonitrile/water (60/40, v/v). Verify that your mobile phase is miscible with the shipping solvent before equilibrating your column for use.

- Flush column for 5 minutes with 100 % acetonitrile at a flow rate of:
  - 0.6 mL/min for analytical dimensions (i.e. – 2.0 mm ID)
  - 1 mL/min for analytical dimensions (i.e. – 3.0 mm ID)
  - 4 mL/min for analytical dimensions (i.e. – 4.6 mm ID).
  - 10 mL/min for semi-prep dimensions (i.e. – 10.0 mm ID)
- Continue conditioning the column with your mobile phase until you get a stable baseline

### **SPECIFICALLY for 150 x 0.1 mm dimension**

Onyx reversed phase 150 x 0.1 mm columns are shipped in methanol/water (80:20).

- Verify that your mobile phase is miscible with the shipping solvent before equilibrating your column for use
- Install your column as described above. Make sure that no air bubbles are in the system
- Equilibrate by passing 10 column volumes of mobile phase at normal flow rate until you achieve a stable baseline

### **SPECIFICALLY for Silica (Si) phase**

Onyx normal phase columns are shipped in n-heptane/dioxane (95/5, v/v). Verify that your mobile phase is miscible with the shipping solvent before equilibrating your column for use. As it can dry out during stocking and shipping thoroughly activate the packing by equilibrating.

- Flush column for 5 minutes with n-heptane/ dioxane (50/50, v/v) at a flow rate of 3 mL/min
- Continue conditioning your column with your mobile phase until you get a stable baseline

## **CLEANING & REGENERATION PROCEDURE**

For cleaning and regeneration of Onyx materials, connect the column with the flow arrow on the label pointing toward the pump for backflushing.

*(Continued on next page)*

## **SPECIFICALLY for reversed phase materials (C18 & C8)**

In most cases, a flushing with 100 % acetonitrile or methanol for 5 minutes (see Table 7) is sufficient. If buffers have been used, first pump 100 % water, then methanol.

<i>TABLE 7</i>	
<b>Column ID (mm)</b>	<b>Flow Rate (mL)</b>
2.0	0.6
3.0	1.0
4.6	3.0
10.0	15.0

If the result is not satisfactory, flush the Onyx reversed phase column (see Table 7) with the following solvents, one after the other, for 5 minutes each in the following order:

1. Water
2. Acetonitrile
3. 2-Propanol
4. Heptane
5. 2-Propanol
6. Acetonitrile
7. Water

## **SPECIFICALLY for 150 x 0.1 mm dimension**

A shift in retention or resolution or unspecific background may indicate contamination of the column.

- Use 95 % acetonitrile for cleaning
- Make sure that your in-column solvent or mobile phase is miscible with the cleaning solvent
- Flush the column with 2 – 4 column volumes of 95 % acetonitrile

## **SPECIFICALLY for Silica (Si) phase**

Flush the Onyx Silica (Si) column (see Table 7) with the following solvents for 5 minutes each in the following order:

- n-heptane
- n-heptane/ dioxane (50/50)
- dioxane
- n-heptane/ dioxane (50/50)

DO NOT use more than 5 % DMSO, 5 % chlorinated hydrocarbons, or solvent mixtures containing more than 50 % THF.

---

## **COLUMN STORAGE**

- Column storage conditions affect the column lifetime
- When storing the column for several days or longer, store the column in 100 % acetonitrile
- Never store column with buffers
- If the mobile phase contained a buffer, flush with 10 column volumes of HPLC grade water to remove any buffers or salts, then with acetonitrile
- Confirm that the column end plugs are firmly in place

## **SPECIFICALLY for 150 x 0.1 mm dimension**

- For prolonged storage, flush the column with a mobile phase with 60 to 80 % acetonitrile or methanol in water
- In case the column has been used with buffer media, flush the column with several column volumes of 60 – 80 % acetonitrile or methanol
- Never store columns for a long time with buffer or acid containing solvents
- When not in use, store the column in the protective shipping box

## **SPECIFICALLY for Silica (Si) phase**

- When storing a normal phase column for several days or longer, store the column in n-heptane/ dioxane (95/5, v/v).

## PART III - SFC (SUPERCRITICAL) FLUID CHROMATOGRAPHY

Phenomenex analytical and Axia 'SFC Approved' columns have been leak tested under SFC conditions at pressure far exceeding what may be expected with normal SFC operation.

---

### RUNNING PARAMETERS



- Backpressure Limitation: 3500 psi
- Flow Rate Limitations: Flow rate is to be controlled so that pressure limit of 3500 psi is not exceeded
- pH Limitations: Dictated by the media packed in the column

---

### EQUILIBRATING COLUMN

SFC column stationary phases have a polar surface and may be shipped under reversed phase or normal phase conditions. Flush all columns with 10-30 column volumes of Methanol/CO<sub>2</sub> as intermediate solvent between CO<sub>2</sub> and column shipping conditions. Be aware of backpressure settings.

Equilibrate column to starting conditions with 10 column volumes of mobile phase.

---

### MOBILE PHASE CONSIDERATIONS

- Use only HPLC grade solvent modifiers
- Use only highest purity chemicals and reagents
- Filter and degas all mobile phases prior to use

---

### CLEANING PROCEDURE

- Under extreme conditions the column can be flushed with 50/50 Acetonitrile/Isopropyl Alcohol followed by 100 % Isopropyl Alcohol. Maintain backpressure below limits.
- Re-Equilibrate column to starting conditions with 10 column volumes of mobile phase

---

### COLUMN STORAGE

- Completely remove all buffers, acids, bases or other mobile phase additives to prevent damage to media
- Flush with at least 10 column volumes of Methanol after the last sample is purified
- Store column with end plugs firmly seated in endfittings to ensure storage solvent does not evaporate

## PART IV - AXIA PACKED PREPARATIVE COLUMN

---

### RUNNING PARAMETERS

- Backpressure Limitation: 3500 psi
- Flow Rate Limitations: Determined by the viscosity of mobile phase; flow rates to be controlled so that backpressure limit of 3500 psi is not exceeded
- pH Limitations: Dictated by the media packed in the column

---

### MOBILE PHASE CONSIDERATIONS

- Use only HPLC grade solvents and water
- Use only highest purity chemicals and reagents
- Filter and degas all mobile phases prior to use

---

### CLEANING PROCEDURE

(for Axia reversed phase columns)

- For achiral applications, see silica-based and twin technology column cleaning procedures.

- Under extreme conditions, the column can be flushed with 10 column volumes of 100 % THF (or IPA) followed by 100 % methylene chloride
- After cleaning, wash with 100 % THF (or IPA) and 50:50 Acetonitrile/Water, prior to equilibrating with the starting mobile phase n.

---

## COLUMN STORAGE

- Completely remove all buffers, acids, bases, or other mobile phase additives to prevent physical damage to the media
- Flush with at least 10 column volumes of 50:50 Acetonitrile/Water after the last sample is purified
- Store with column end plugs placed back in the end-fittings to ensure that the packing media does not dry out

*For additional information, consult the Care and Use of Axia Packed Preparative HPLC Columns, included with each Axia column purchased.*

## PART V - LUX CHIRAL COLUMNS

### RUNNING PARAMETERS

#### OPERATING BACKPRESSURE

The mobile phase flow rate should be set such that the column backpressure stays below 300 bar (4300 psi). This maximum backpressure should not be exceeded for long periods of time.

#### OPERATING TEMPERATURES

With standard mobile phases (such as alkane/alcohol) the column can be used in the temperature range 0-50 °C.

---

### MOBILE PHASE CONSIDERATIONS

#### MOBILE PHASE COMPATIBILITY

Lux columns can be used with normal phase (alkane/alcohol), reversed phase (aqueous methanol, aqueous acetonitrile or appropriate buffer/methanol or buffer/acetonitrile mixtures), as well as with pure polar organic solvents (low molecular weight alcohols, acetonitrile or their mixtures).

#### SOLVENT SWITCHING

An appropriate column washing procedure must be applied when changing from one mobile phase to another. The miscibility of the different mobile phase components must be carefully considered for this wash.

To safely transfer a column from normal phase to polar organic or reversed phase conditions, flush the column with methanol/ethanol 9:1 (V/V) as transition solvent at a flow rate of 0.5 mL/min. Flush the column with at least ten column volumes (i.e. 25 mL for a 250 x 4.6 mm ID column or 15 mL for a 150 x 4.6 mm ID column) to completely remove the initial mobile phase. When the column has been flushed, equilibrate the column with at least ten column volumes of the polar organic or reversed phase solvent mixture to condition the column. In addition, when the buffer salt additive of the reversed phase mobile phase is insoluble in methanol/ethanol, flush the column briefly with water before switching to the buffered mobile phase. When the column has been flushed equilibrate the column with at least ten column volumes of the reversed phase solvent mixture.

To safely transfer a column from polar organic to normal phase conditions flush the column with at least ten column volumes of methanol/ethanol 9:1 (V/V) as transition solvent at a flow rate of 0.5 mL/min. When the column has been flushed with methanol/ethanol equilibrate the column with at least ten column volume of the normal phase solvent mixture to condition the column. We do not recommend switching from reversed phase mode back to normal phase mode.

## USE OF MOBILE PHASE MODIFIERS

For basic samples or acidic chiral compounds, it may be necessary to use an appropriate mobile phase modifier in order to achieve chiral resolution and to insure proper peak shapes. Diethylamine, ethanolamine and butyl amine in the concentration range 0.1-0.5 % can be used with basic analytes, while trifluoroacetic or acetic acid (0.1-0.5 %; typically 0.1-0.2 %) with acidic analytes. Mixtures of basic and acidic mobile phase additives are acceptable (e.g. diethylamine acetate or trifluoroacetate). Lux columns will deliver consistent results when operated with mobile phases containing additives at the concentration levels specified above. However, limited decrease in column efficiency may occur when a column is used in combination with these additives. Therefore, we advise to dedicate columns to mobile phases containing basic additives.

## MOBILE PHASE RESTRICTIONS

Lux chiral stationary phases are prepared by coating silica with various polysaccharide derivatives. Therefore, any solvent dissolving the polysaccharide derivative (such as tetrahydrofuran, acetone, chlorinated hydrocarbons, ethylacetate, dimethylsulfoxide, dimethylformamide, N-methylformamide, toluene, methylethyl ketone and methyl tert-butyl ether, etc. must be avoided even in trace amounts (e.g. even as sample solvent).

---

## EXTENDING LIFETIME AND RECONDITIONING

Phenomenex recommends the use of SecurityGuard™ guard cartridges to extend the lifetime of your column, especially with samples extracted from complex matrixes. Ideally, samples must be completely dissolved in the mobile phase or filtered through a syringe filter of approximately 0.45 µm porosity.

To regenerate or remove potential contaminant after extended use of your Lux column, we recommend flushing the column with methanol for polar organic and reversed phase mode or with ethanol for normal phase mode for 2-3 hours at the appropriate flow rate.

---

## COLUMN STORAGE

- Column storage for a longer period of time is recommended in n-hexane/2-propanol (9:1, v/v).
- Columns used in reversed phase conditions should be first flushed with water (whenever a buffer salt was used as RP mobile phase additive) and then with methanol (or with methanol only when no salt was used). The column can be stored in methanol

# PART VI - CHIREX CHIRAL COLUMNS

## RUNNING PARAMETERS

- Temperature must not exceed 50 °C
- Column pressure must not exceed 3000 psi
- Maintain flow rate between 0.5-2.0 mL/min for 4.6 mm ID columns

---

## MOBILE PHASE CONSIDERATIONS

- Dedicate column to reversed or normal phase solvents
- pH range: 2.5 to 7.5
- Use only HPLC grade solvents
- Use only highest purity chemicals and reagents

- Filter and degas all mobile phases prior to use
- Make sure solvents are miscible (see pp. 22-23)

*Most CHIREX Chiral columns use a Type I or brush type Chiral stationary phase (CSP). Normal phase systems usually provide better selectivity than reversed phase systems. SEE COLUMN INSERT FOR FURTHER INFORMATION ABOUT SPECIFIC CHIREX COLUMNS (included with each column)*

---

## SAMPLE CONSIDERATIONS

Always prefilter samples with Phenex™ syringe filters to avoid particulate contaminants which may clog the column. Use of a guard column is highly recommended to prolong the life of your analytical to preparative column.

## PART VII - BIOSEP & YARRA SEC COLUMNS

### RUNNING PARAMETERS

- Maximum flow rate: 1.5 mL/min
- Column pressure must not exceed 1500 psi for BioSep and 3000 psi for Yarra SEC-2000 or -3000
- Column pressure for Yarra SEC-4000 must not exceed 1750 psi
- Maximum temperature: 50 °C

---

### MOBILE PHASE CONSIDERATIONS

- pH range: 2.5 - 7.5
- Maximum organic modifier: Up to 100 % CH<sub>3</sub>CN, 10 % DMSO or 500 mM β-mercaptoethanol
- Maximum salt concentration: 1 M
- Filter and degas all mobile phases prior to use

---

### SAMPLE CONSIDERATIONS

Always prefilter samples with Phenex™ syringe filters and add SecurityGuard™ guard cartridge system to maximize column lifetime.

---

### CLEANING PROCEDURE

- General protein removal: wash with 30 mL of 0.1 M NaH<sub>2</sub>PO<sub>4</sub>, pH 3.0
- Hydrophobic protein removal: use Acetonitrile gradient
- Strongly adsorbed proteins: wash with 30 mL of 6 M guanidine thiocyanate or 10 % DMSO

---

### COLUMN STORAGE

- Overnight storage: run mobile phase at 0.2 mL/min
- Prolonged storage: use 0.05 % sodium azide in water or 20 % methanol in water
- Yarra SEC-4000 columns should only be stored in 20 % methanol. Be careful not to exceed 1750 psi when exchanging between buffer and storage conditions.

## PART VIII - REZEX POLYMER-BASED COLUMN

### RUNNING PARAMETERS

- Columns should be run at elevated temperatures
- (60-85 °C) except Rezex ROA and RHM for most applications (Rezex ROA and RHM ~ 40 °C)
- Column pressure for 8 % cross-linked material must not exceed 1,000 psi; must not exceed 300 psi for 4 % cross-linked material



- Clean and reverse flush column regularly with HPLC grade water
- To increase column lifetime, use a Rezex guard column or SecurityGuard™ cartridge system (See p.16-19 and p.33)

**Important: Never exceed maximum pressure limitations. This will cause irreversible damage to the column.**

## MOBILE PHASE CONSIDERATIONS

- Filter and degas all mobile phases prior to use
- Replace mobile phase frequently to avoid microbial contamination
- Do not exceed 5 % Methanol, IPA, EtOH
- Do not exceed 30 % Acetonitrile or other organic
- Store columns in HPLC grade water

Rezex utilizes a sulfonated polystyrene resin which is very rugged and resistant to chemical attack. However, the material **is pressure sensitive and must be cared for properly.**

### START UP (IMPORTANT!)

Turn on column heating unit to 60 - 85 °C and start the mobile phase at 0.1 mL/min. Make sure the pressure remains below 400 psi for 8 % cross-linked material; below 200 psi for 4 % cross-linked material. As the temperature reaches working condition, increase flow rate to the specified level. (See Rezex Operating Parameters)

Column may exhibit brown/gray liquid upon startup after storage. This is perfectly normal and will dissipate after a few minutes. It will not adversely affect your LC system.

### SHUT DOWN

- Reduce flow rate slowly

Overnight: Lower flow rate to 0.1 mL/min. Leave system on and continue heating.

Long Term: Store columns in 100 % water. Turn off pump and allow the system to cool. Replace the end plugs and tightly cap the column.

## SAMPLE CONSIDERATIONS

Always prefilter samples with Phenex™ syringe filters to avoid particulate contaminants which may clog the column. Use of a guard column is highly recommended to prolong the life of your analytical column.

## CLEANING PROCEDURE

Before utilizing any cleaning procedure outlined in the Tables on pages 17 and 19, first try to clean your Rezex column as follows:

Remove the guard column and reverse the direction of flow on the analytical column. Run 100 % HPLC grade water through the column as follows:

TABLE 8		
Rezex Column	Flow (mL/min)	Temp. (°C)
RPM, RCM, RHM	0.4	85
RCU	0.2	85
RSO and RNO	0.1	75
RNM and RAM	0.4	75
ROA	0.4	85

Run the column under these conditions for a minimum of 12 hours. After completing the cleaning procedure, return the column to the original direction of flow and equilibrate for analysis.

If this procedure is not effective in cleaning the column, proceed to the specified procedures outlined in Tables 9 and 10.

# PART VIII - REZEX POLYMER-BASED COLUMNS (cont'd)

## SPECIFICATIONS AND OPERATING PARAMETERS

<b>Table 9</b>	<b>RCM Monosaccharide</b>	<b>RSO Oligosaccharide</b>	<b>RNO Oligosaccharide</b>	<b>RNM Carbohydrate</b>	<b>RAM Carbohydrate</b>
<b>Part Number</b>	00H-0130-K0	00P-0133-N0	00P-0137-N0	00H-0136-K0	00H-0131-K0
<b>Ionic Form</b>	Calcium	Silver	Sodium	Sodium	Silver
<b>Standard Dimensions</b>	300 x 7.8 mm	200 x 10 mm	200 x 10 mm	300 x 7.8 mm	300 x 7.8 mm
<b>Matrix</b>		Sulfonated Styrene Divinyl Benzene			
<b>Cross Linking</b>	8 %	4 %	4 %	8 %	8 %
<b>Particle Size (µm)</b>	8	12	12	8	8
<b>Min. Efficiency (p/m) based on last peak</b>	35,000	N/A	N/A	30,000	35,000
<b>Typical Pressure (psi @ Testing Flow Rate)</b>	260	115	130	170	285
<b>Max. Pressure (psi @ Max Flow Rate)</b>	1,000	300	300	1,000	1,000
<b>Max. Flow Rate (mL/min)*</b>	1.0	0.3	0.3	1.0	1.0
<b>Max. Temperature (°C)</b>	85	85	85	85	85
<b>Typical Mobile Phase</b>	Water	Water	Water	Water	Water
<b>pH Range</b>	Neutral	Neutral	Neutral	Neutral	Neutral
<b>Guard Column Part No.</b>	03B-0130-K0	03R-0133-N0	03R-0137-N0	03B-0136-K0	03B-0131-K0

\* Make sure the maximum pressure is not exceeded

# COLUMNS FOR CARBOHYDRATE AND ORGANIC ACID ANALYSIS

**Table 9 (continued)**

<b>Cleaning, Regeneration and Storage</b>		<b>RCM Monosaccharide</b>	<b>RSO Oligosaccharide</b>	<b>RNO Oligosaccharide</b>	<b>RNM Carbohydrate</b>	<b>RAM Carbohydrate</b>
Organic Modifiers (Max)			5 % Methanol, IPA, EtOH			
Inorganic Modifiers (Max)		5 % CaSO <sub>4</sub> , Ca(NO <sub>3</sub> ) <sub>2</sub> , CaCl <sub>2</sub>	5 % Silver Nitrate	5 % Sodium Salts	5 % Sodium Salts	2 % Silver Nitrate
Avoid		Acids, Bases, Non-Calcium Salts or Metal Ions, >30 % Acetonitrile	Acids, Bases, Non-Silver Salts/Metal Ions, >30 % Acetonitrile	Acids, Bases, Non-Sodium Salts/Metal Ions, >30 % Acetonitrile	Acids, Bases, Non-Sodium Salts/Metal Ions, >30 % Acetonitrile	Acids, Bases, Non-Silver Salts/Metal Ions, >30 % Acetonitrile
Cleaning Solvent		100 % Water	100 % Water	100 % Water	100 % Water	100 % Water
Flow Rate (mL/min)		0.4	0.1	0.4	0.4	0.4
Temperature (°C)		85	85	85	85	85
Duration (hrs)		12	12	12	12	12
Regeneration Solvent		0.1 M Ca(NO <sub>3</sub> ) <sub>2</sub>	0.1 M AgNO <sub>3</sub>	0.1 M NaNO <sub>3</sub>	0.1 M NaNO <sub>3</sub>	0.1 M AgNO <sub>3</sub>
Flow Rate (mL/min)		0.2	0.1	0.1	0.2	0.2
Temperature (°C)		85	85	85	85	85
Duration (hrs)		4-16	4-16	4-16	4-16	4-16
Ship/Storage Solvent		Water	Water	Water	Water	Water

# PART VIII - REZEX POLYMER-BASED COLUMNS (cont'd)

## SPECIFICATIONS AND OPERATING PARAMETERS

<b>Table 10</b>	<b>RPM Monosaccharide</b>	<b>RHM Monosaccharide</b>	<b>ROA Organic Acid</b>	<b>RFQ Fast Acid</b>	<b>RCU Sugar Alcohols</b>
<b>Part Number</b>	00H-0135-K0	00H-0132-K0	00H-0138-K0	00D-0223-K0	00G-0130-D0
<b>Ionic Form</b>	Lead	Hydrogen	Hydrogen	Hydrogen	Calcium
<b>Standard Dimensions</b>	300 x 7.8 mm	300 x 7.8 mm	300 x 7.8 mm	100 x 7.8 mm	250 x 4.0 mm
<b>Matrix</b>	Sulfonated Styrene Divinyl Benzene				
<b>Cross Linking</b>	8 %	8 %	8 %	8 %	8 %
<b>Particle Size (µm)</b>	8	8	8	8	8
<b>Min. Efficiency (p/m) (based on last peak)</b>	35,000	35,000	50,000 (Acetic Acid)	30,000	12,000
<b>Typical Pressure (psi @ Testing Flow Rate)</b>	190	275	580	365	90
<b>Max. Pressure (psi @ Max Flow Rate)</b>	1,000	1,000	1,000	1,000	1,000
<b>Max. Flow Rate (mL/min)*</b>	1.0	1.0	1.0	1.0	0.5
<b>Max. Temperature (°C)</b>	85	85	85	85	85
<b>Typical Mobile Phase</b>	Water	Water	0.005N H <sub>2</sub> SO <sub>4</sub>	0.005N H <sub>2</sub> SO <sub>4</sub>	Water
<b>pH Range</b>	Neutral	1-8	1-8	1-8	Neutral
<b>Guard Column Part No.</b>	03B-0135-K0	03B-0132-K0	03B-0138-K0	03B-0223-K0	03A-0130-D0

\* Make sure the maximum pressure is not exceeded

# COLUMNS FOR CARBOHYDRATE AND ORGANIC ACID ANALYSIS

**Table 10 (continued)**

<b>Cleaning, Regeneration and Storage</b>	<b>RPM Monosaccharide</b>	<b>RHM Monosaccharide</b>	<b>ROA Organic Acid</b>	<b>RFQ Fast Acid</b>	<b>RCU Sugar Alcohols</b>
Organic Modifiers (Max)	5 % Lead Nitrate	5 % HNO <sub>3</sub> , H <sub>3</sub> PO <sub>4</sub>	5 % Methanol, IPA, EtOH	5 % HNO <sub>3</sub> , H <sub>3</sub> PO <sub>4</sub>	5 % CaSO <sub>4</sub> , Ca(NO <sub>3</sub> ) <sub>2</sub> , CaCl <sub>2</sub>
Inorganic Modifiers (Max)	Acids, Bases, Non-Lead Salts/Metal Ions, >30 % Acetonitrile	Acids, Bases, Salts, Metal Ions, pH > 3, >30 % Acetonitrile	Acids, Bases, Salts, Metal Ions, pH > 3, >30 % Acetonitrile	Acids, Bases, Salts, Metal Ions, pH > 3, >30 % Acetonitrile	Acids, Bases, Non-Calcium Salts or Metal Ions, >30 % Acetonitrile
Cleaning Solvent	100 % Water	100 % Water	100 % Water	100 % Water	100 % Water
Flow Rate(mL/min)	0.4	0.4	0.4	0.4	0.1
Temperature (°C)	85	85	85	85	85
Duration (hrs)	12	12	12	12	12
Regeneration Solvent	0.1 M Pb(NO <sub>3</sub> ) <sub>2</sub>	0.025 M H <sub>2</sub> SO <sub>4</sub>	0.025 M H <sub>2</sub> SO <sub>4</sub>	0.025 M H <sub>2</sub> SO <sub>4</sub>	0.1 M Ca (NO <sub>3</sub> ) <sub>2</sub>
Flow Rate (mL/min)	0.2	0.2	0.2	0.2	0.1
Temperature (°C)	85	85	85	85	85
Duration (hrs)	4-16	4-16	4-16	4-16	4-16
Ship/Storage Solvent	Water	Water	0.005 N H <sub>2</sub> SO <sub>4</sub>	0.005 N H <sub>2</sub> SO <sub>4</sub>	Water

See pp. 14-15 for general care and usage of Rezex columns.

# PART IX - POLYSEP-GFC-P COLUMNS

## RUNNING PARAMETERS

- Column pressure must not exceed 1000 psi
- Do not exceed 60 °C

---

## MOBILE PHASE CONSIDERATIONS

- pH range: 3 - 12
- Maximum salt concentration: 0.5 M
- Organic Modifier capacity:

	<i>POLYSEP PHASE</i>						
	<b>1000</b>	<b>2000</b>	<b>3000</b>	<b>4000</b>	<b>5000</b>	<b>6000</b>	<b>Linear</b>
<b>Methanol</b>	20 %	95 %	70 %	70 %	70 %	70 %	70 %
<b>Acetonitrile</b>	20 %	70 %	70 %	70 %	70 %	70 %	70 %

---

## CLEANING PROCEDURE

0.5 % SDS or 6 M guanidine thiocyanate. All PolySep columns except for PolySep 1000 may also be cleaned with 50 % acetonitrile. Make sure not to exceed a maximum pressure of 650 psi when cleaning.

---

## COLUMN STORAGE

- Overnight storage: run water at low flow rate (0.2 mL/min or less)
- Prolonged storage: store in 0.05 % sodium azide in water or 10 % methanol in water

---

## SAMPLE CONSIDERATIONS

Always prefilter samples with Phenex™ syringe filters to avoid particulate contaminants which may clog the column. Use of a guard column is highly recommended to prolong the life of your analytical to preparative column.

# PART X - PHENOGEL GPC COLUMNS

## SPECIFICATIONS

<b>Matrix:</b>	Styrene-Divinyl Benzene Copolymer
<b>Particle Size:</b>	5, 10 $\mu\text{m}$
<b>Porosities:</b>	50 $\text{\AA}$ to $10^6$ $\text{\AA}$ , and mixed beds
<b>Maximum Pressure***:</b>	1500 psi
<b>Maximum Temperature:</b>	140 $^{\circ}\text{C}$
<b>Minimum Efficiency*:</b>	<b>5 <math>\mu\text{m}</math>:</b> 45,000 P/m** <b>10 <math>\mu\text{m}</math>:</b> 35,000 P/m**
<b>Typical Flow Rates:</b>	<b>4.6 mm ID:</b> 0.35 mL/min <b>7.8 mm ID:</b> 1.0 mL/min <b>21.2 mm ID:</b> 7.0 mL/min
<b>End Fittings:</b>	Valco Compatible

\*Tested in THF      \*\* For 300 x 7.8 mm ID columns

\*\*\* At testing flow rates in THF, other solvents will vary

## SAMPLE CONSIDERATIONS

Always prefilter samples with Phenex™ syringe filters to avoid particulate contaminants which may clog the GPC column. Use of a GUARD COLUMN is highly recommended to prolong the life of your analytical or preparative column. For optimal results, use the chart below to determine sample concentrations and injection volumes.

TABLE 11

Molecular Weight	Concentration (w/v)	Max Injection Volume
< 50 K	0.5 %	100 $\mu\text{L}$
50-600 K	0.25 %	100 $\mu\text{L}$
600-3000 K	0.05 %	100 $\mu\text{L}$
>3000 K	0.01 %	20 $\mu\text{L}$

Continued on p. 24

# SOLVENT MISCIBILITY TABLE

TABLE 12

## Solvent Polarity Chart

Relative Polarity	Compound Formula	Group	Representative Solvent Compounds
<b>NONPOLAR</b>	R - H	Alkanes	Petroleum ethers, ligroin, hexanes
	Ar - H	Aromatics	Toluene, benzene
	R - O - R	Ethers	Diethyl ether
	R - X	Alkyl halides	Tetrachloromethane, chloroform
	R - COOR	Esters	Ethyl acetate
	R - CO - R	Aldehydes and ketones	Acetone, methyl ethyl ketone (MEK)
	R - NH <sub>2</sub>	Amines	Pyridine, triethylamine
	R - OH	Alcohols	Methanol, ethanol, isopropanol, butanol
	R - COHN <sub>2</sub>	Amides	Dimethylformamide
	R - COOH	Carboxylic acids	Ethanoic acid
<b>POLAR</b>	H - OH	Water	Water

Increasing Polarity

Solvent	Polarity Index	Refractive Index @ 20 °C	UV (nm) Cutoff @ 1 AU	Boiling Point (°C)	Viscosity (cPoise)	Solubility in Water (% w/w)
Acetic Acid	6.2	1.372	230	118	1.26	100
Acetone	5.1	1.359	330	56	0.32	100
Acetonitrile	5.8	1.344	190	82	0.37	100
Benzene	2.7	1.501	280	80	0.65	0.18
Butyl Acetate	4.0	1.394	254	125	0.73	0.43
n-Butanol	3.9	1.399	215	118	2.98	7.81
Carbon tetrachloride	1.6	1.466	263	77	0.97	0.08
Chloroform	4.1	1.446	245	61	0.57	0.815
Cyclohexane	0.2	1.426	200	81	1.00	0.01
1,2-Dichloroethane <sup>1</sup>	3.5	1.444	225	84	0.79	0.81
Dichloromethane <sup>2</sup>	3.1	1.424	235	41	0.44	1.6
Dimethylformamide	6.4	1.431	268	155	0.92	100
Dimethyl sulfoxide <sup>3</sup>	7.2	1.478	268	189	2.00	100
Dioxane	4.8	1.422	215	101	1.54	100
Ethyl Acetate	4.4	1.372	260	77	0.45	8.7
Ethanol	5.2	1.360	210	78	1.20	100
di-Ethyl Ether	2.8	1.353	220	35	0.32	6.89
Heptane	0.0	1.387	200	98	0.39	0.0003





## COLUMN STORAGE

Solvents such as THF (stabilized THF only), Chloroform, Methylene Chloride, and Toluene are commonly used for column storage. Be sure to follow solvent switching instructions (see below) if using solvents other than THF. Storage solvents that remain liquefied at ambient temperatures and are not oxidizing can be used for storage.

**BE SURE THAT ANY COLUMN THAT IS NOT USED IS CAPPED TIGHTLY WITH END-PLUGS TO AVOID EVAPORATION OF SOLVENTS FROM COLUMN. COLUMN DESICCATION IS THE MOST COMMON SOURCE OF COLUMN FAILURE.**

## SOLVENT SWITCHING CONSIDERATIONS FOR NON-AQUEOUS GPC COLUMNS

Phenogel columns are rugged and exhibit wide solvent compatibility. Different solvents, however, produce different swell characteristics (Table 13). Improper solvent switches can result in a void. For this reason, we recommend that you dedicate columns to specific solvents.

If you need to switch solvents, it is VERY IMPORTANT to take the following into consideration:

1. Reduce flow rate to 0.2 mL/min.
2. Backpressure must NEVER exceed 650 psi.
3. Always check solvent miscibility in a beaker or follow the solvent miscibility table on page 22-23 before proceeding with ANY solvent switch.
4. Compare the swell characteristics of solvent 1 (old solvent) to solvent 2 (new solvent) and use the following guidelines:
  - If solvent 1 and solvent 2 belong to the same swell category (Table 13), check the solvent miscibility and proceed with the switch.
  - If solvent 1 and solvent 2 belong to successive swell categories as indicated by the arrows on Table 13, check the miscibility and proceed with the switch.
  - If solvent 1 and solvent 2 DO NOT belong to the same OR successive swell categories, switch to an intermediate solvent FIRST, as indicated by the arrows on Table 13.

Table 13		
A – 70-80 % SWELL		
Toluene	TCB	
THF	p-Xylene	
Chloroform	Cyclopentane	
Benzene	Pyridine	
o-Dichlorobenzene		

B – 60 % SWELL	
Ethyl Ether	Cycloheptane
Methylene Chloride	Ethyl Acetate
Methyl Ethyl Ketone	

C – 50 % SWELL	
Acetone	Dimethyl Acetamide
m-Cresol	N-Methyl Pyrrolidone
o-Chlorophenol	DMSO
Dimethyl Formamide	Dioxane

D – 30 % SWELL		
Acetonitrile	Freon 113	n-Butyl Alcohol
Cyclohexane	HFIP	Methanol
n-Hexane	Iso-Propanol	

# SOLVENT COMPATIBILITY CHART FOR PHENOGEL GPC COLUMNS

Mobile Phase Solvent	PHENOGEL PORE SIZE								Suggested Operating Temperature
	50 Å	100 Å	500 Å	10 <sup>3</sup> Å	10 <sup>4</sup> Å	10 <sup>5</sup> Å	10 <sup>6</sup> Å	Linear & Mixed	
Acetone	Y	Y	Y	Y	Y	Y	Y	Y	
Benzene	Y	Y	Y	Y	Y	Y	Y	Y	
Carbon Tetrachloride	Y	Y	Y	Y	Y	Y	Y	Y	
Chloroform	Y	Y	Y	Y	Y	Y	Y	Y	
30 % HFIP/chloroform	Y	Y	Y	Y	Y	Y	Y	Y	
Diethyl Ether	Y	Y	Y	Y	Y	Y	Y	Y	
Dimethylacetamide (DMAC)	Y*	Y	Y	Y	Y	Y	Y	Y	60 °C
Dimethylformamide (DMF)	Y*	Y	Y	Y	Y	Y	Y	Y	60 °C
Dioxane	Y	Y	Y	Y	Y	Y	Y	Y	
DMSO	Y*	Y	Y	Y	Y	Y	Y	Y	60 °C
Ethyl Acetate	Y	Y	Y	Y	Y	Y	Y	Y	
Hexafluoroisopropanol (HFIP)	Y	Y	Y	Y	Y	Y	Y	Y	

\*Not recommended on 5 µm 50 Å columns.

N = Not Compatible Y = Compatible

Table 14

Mobile Phase Solvent	PHENOGEL PORE SIZE								Suggested Operating Temperature
	50 Å	100 Å	500 Å	10 <sup>3</sup> Å	10 <sup>4</sup> Å	10 <sup>5</sup> Å	10 <sup>6</sup> Å	Linear & Mixed	
Hexane	Y	Y	Y	Y	Y	Y	Y	Y	
m-Cresol	Y*	Y	Y	Y	Y	Y	Y	Y	100 °C
Methyl Ethyl Ketone	Y	Y	Y	Y	Y	Y	Y	Y	
Methylene Chloride	Y	Y	Y	Y	Y	Y	Y	Y	
o-Chlorophenol	Y*	Y	Y	Y	Y	Y	Y	Y	100 °C
o-Dichlorobenzene	Y*	Y	Y	Y	Y	Y	Y	Y	135 °C
Quinolin	Y*	Y	Y	Y	Y	Y	Y	Y	60 °C
Tetrahydrofuran	Y	Y	Y	Y	Y	Y	Y	Y	
Toluene	Y	Y	Y	Y	Y	Y	Y	Y	
Trichlorobenzene	Y*	Y	Y	Y	Y	Y	Y	Y	135 °C
Water	N	N	N	N	N	N	N	N	
Xylene	Y	Y	Y	Y	Y	Y	Y	Y	

## PART XI - POLYMERX RP COLUMNS

### SPECIFICATIONS

- Matrix: Polystyrene Divinylbenzene (PSDVB)
- Particle Size: 3, 5, 7, 10  $\mu\text{m}$
- Pore Size: 100  $\text{\AA}$

### RUNNING PARAMETERS

- Maximum temperature: 60  $^{\circ}\text{C}$
- Maximum pressure: 2500 psi

### MOBILE PHASE CONSIDERATIONS

- pH range: 0 - 14
- Avoid buffer strength  $> 0.5 \text{ N}$

### CLEANING PROCEDURE

- 100 % Water to 100 % Acetonitrile. Repeat 3 times.

### COLUMN STORAGE

- 75:25 Acetonitrile / Water

## PART XII - HPLC COLUMN PROTECTION & PERFORMANCE TESTING

- Maximize the life of your valuable HPLC Column
- Reduce system wear and tear
- Save time and money

### PHENEX™ SYRINGE FILTERS

- Increase column lifetime (save money!)
- Ensure more accurate, consistent results
- Eliminate damaging microparticulates

Particulates can damage expensive equipment, valves, columns and pumps. They can also lead to erratic analytical results. Pre-filtering samples prior to analysis is critical in preventing column and frit blockage, undue wear on valve seals, and abnormally high operating pressures.



TABLE 15

Sample or Mobile Phase Volume (mL)	Filter Membrane (diameter, mm)	Format
$\leq 2$	4	Syringe filter
2 to 10	15	Syringe filter
10 to 100	25-28	Syringe filter
$> 100$	47	Membrane disk
$> 1000$	90	Membrane disk

### MEMBRANE FILTERS ORDER LIST GUIDE

#### REGENERATED CELLULOSE (RC)

As a universal hydrophilic membrane, RC is widely used in chromatography for the clarification of aqueous samples and solvents. Due to its ultra-low binding capabilities, RC membranes are an excellent choice for proteins, peptides and other biomolecules.

#### POLYTETRAFLUOROETHYLENE (PTFE, TEFLON®)

PTFE is an inherently hydrophobic membrane, excellent for filtration of organic-based, highly acidic or basic samples and solvents. Widely used in chromatography, it is especially well suited for the clarification of non-aqueous samples. Although this membrane is hydrophobic, it can be made hydrophilic by

wetting the membrane with alcohol and then flushing with deionized water.

## POLYETHERSULFONE (PES)

Polyethersulfone, a hydrophilic membrane with fast flow, high-throughput characteristics, with ultra-low protein binding. It is ideally suited for use in life sciences applications. The PES membrane offers better chemical resistance than cellulose acetate. Recommended for filtering critical biological samples, tissue culture media, additives, and buffers.

## NYLON (NY)

Nylon has inherent hydrophilic characteristics and works well for filtration of many aqueous and mixed-organic samples. Nylon exhibits a high non-specific affinity for proteins. Phenomenex recommends Phenex-RC (Regenerated Cellulose) filters for applications requiring low non-specific adsorption of proteins.

## CELLULOSE ACETATE (CA)

Cellulose Acetate membranes exhibit ultra-low protein binding and are broadly used in the filtration of biological samples. In combination with a glass pre-filter (Phenex-GF/CA), this membrane is excellent for filtration of tissue culture media, general biological sample filtration and clarification.

## GLASS FIBER (GF)

Glass Fiber filters are made of inert borosilicate glass and have a nominal 1.2  $\mu\text{m}$  pore size. They are commonly used with highly viscous samples or samples containing high concentrations of particulate matter (e.g., food analysis, biological samples, soil samples, fermentation broth samples, removal of yeasts, molds, etc.). Glass Fiber filters can be used alone or in conjunction with other Phenex filter membranes such as the 0.45  $\mu\text{m}$  pore Phenex-RC filter to reduce clogging of the membrane and optimize flow.

## ORDERING INFORMATION

Part No.	Pore Size ( $\mu\text{m}$ )	Phenex Membrane	Housing
<b>4 mm Diameter (500/pk)</b>			
AF0-3103-52	0.45	RC	PP
AF0-3102-52	0.45	PTFE <sup>6</sup>	PP
AF3-3107-52	0.45	NY	PP
AF0-3203-52	0.20	RC	PP
AF0-3202-52	0.20	PTFE <sup>6</sup>	PP
AF3-3207-52	0.20	NY	PP
<b>15 mm Diameter (500/pk)</b>			
AF0-2103-52	0.45	RC	PP
AF0-2102-52	0.45	PTFE <sup>6</sup>	PP
AF0-2107-52	0.45	NY	PP
AF0-2203-52	0.20	RC	PP
AF0-2202-52	0.20	PTFE <sup>6</sup>	PP
AF0-2207-52	0.20	NY	PP
<b>25–28 mm Diameter (500/pk)</b>			
AF0-8103-52 <sup>5</sup>	0.45	RC	PP
AF0-8108-52 <sup>7</sup>	0.45	PES <sup>3</sup>	PP
AF0-1102-52	0.45	PTFE <sup>6</sup>	PP
AF0-1107-52	0.45	NY	PP
AF0-8B09-52 <sup>7</sup>	0.45	GF/CA <sup>2,3,4</sup>	MBS
AF0-8203-52 <sup>5</sup>	0.20	RC	PP
AF0-8208-52 <sup>7</sup>	0.20	PES <sup>3</sup>	PP
AF0-1202-52	0.20	PTFE <sup>6</sup>	PP
AF0-1207-52	0.20	NY	PP
AF0-8A09-52 <sup>7</sup>	0.20	GF/CA <sup>2,3,4,7</sup>	MBS
AF0-8515-52 <sup>7</sup>	1.20	GF <sup>2,3</sup>	MBS

Housing is made of medical-grade polypropylene (PP), unless otherwise indicated. Above syringe filters are non-sterile.

- Additional membrane types available.
- Glass fiber filters are 28 mm diameter and made of borosilicate. They will remove 90 % of all particles >1.2  $\mu\text{m}$ .
- Housing material is methacrylate butadiene styrene (MBS) polymerisate. Also known as Cryolite®.
- Cellulose acetate is surfactant-free.
- 26 mm diameter.
- Hydrophobic membrane. Can be made hydrophilic by pre-wetting with IPA.
- 28 mm diameter.

---

## PHENEX™ DISPOSABLE CENTRIFUGAL FILTER UNITS

- Convenient filtration of multiple HPLC and GC samples
- High recovery for small samples
- Nylon, Cellulose Acetate, and PTFE (Teflon®) membrane materials



Centrifugal force drives the sample through the filter quickly without effort on the part of the chemist. No cleaning of syringes is required between samples. The receiver tube serves as a container for the filtered sample and can be retained as long as desired.

---

## ORDERING INFORMATION

Part No.	Pore Size (µm)	Volumes (mL) Sample/Receiver	Membrane Non-Sterile	Unit
AF0-0438	0.2	2.0 / 5.0	Nylon	25/pk
AF0-0439	0.45	2.0 / 5.0	Nylon	25/pk
AF0-0440	0.2	2.0 / 5.0	PTFE	25/pk
AF0-0441	0.45	2.0 / 5.0	PTFE	25/pk
AF0-8353	0.2	2.0 / 5.0	CA	25/pk
AF0-8354	0.45	2.0 / 5.0	CA	25/pk

*Above centrifugal filters are non-sterile.*

# GUARD CARTRIDGE SYSTEM



SecurityGuard provides a great balance of convenience, column protection capability and value. If you've ever used another guard cartridge system or conventional guard column, you will be pleasantly surprised when you see how practical and effective SecurityGuard really is. This highly advanced, patented design offers several unique features up to now not available.

## CONVENIENCE



Knowing when to replace your guard is no longer a mystery! SecurityGuard's direct-view feature lets you inspect the packing material for visual contaminants and indicates when it's time to replace the cartridge. No other guard cartridge has this convenient feature.

## EXTRA PROTECTION



SecurityGuard offers the option of stacking two cartridges in the same holder, using the simple stacking ring provided. Extra length provides extra protection. When the first cartridge becomes exhausted, contaminants are retained by the second cartridge.

## VERSATILITY



One direct-connect holder conveniently finger-tightens into virtually any brand of HPLC column worldwide. How can one holder be direct-connect and universal at the same time when end-fittings have different depths? Answer- the length of the stainless steel nib at the end of the holder automatically adjusts to the precise depth of a column's endfitting. SecurityGuard's fingertight connection will withstand pressures up to 3500 psi (241 bar) and it features a completely inert and biocompatible flowpath.

## ACCURACY

The cartridges can be used with virtually any matching phase of virtually any brand of column without affecting efficiency, retention time or backpressure. There are 37 different phases to choose from, including cartridges for general purpose, pharmaceutical, protein and polypeptide, aqueous size exclusion, chiral, carbohydrate and organic acid applications. SecurityGuard phases can be used with columns containing 3, 3.5, 4, 5, 10, 15 µm or larger diameter particle sizes.

## SECURITYGUARD ORDERING INFORMATION

### Analytical Holder Assembly Kit

Part No.	Description	Unit
<b>KJ0-4282</b>	Guard Cartridge Kit	ea




*Kit includes: 1 Cartridge Holder, 3 PEEK Ferrules, 2 Stacking Rings, 2 PEEK Fingertight Male Nuts, 2 Wrenches*

*(Continued on next page)*

# SECURITYGUARD ORDERING INFORMATION (CONTINUED)

## Semi-Preparative and Preparative Holder for 10.0, 21.2 and 30.0 mm ID cartridges

Part No.	Description	Unit
<b>AJO-7220</b>	Holder for 10.0 mm ID cartridges	ea
<b>AJO-8223</b>	Holder for 21.2 mm ID cartridges	ea
<b>AJO-8277</b>	Holder for 30.0mm ID cartridges	ea



## Cartridges - General Purpose / Pharmaceutical

Part No.	Material Description	pH Stability	Dimensions L x ID(mm)	Unit
<b>AJO-4286</b>	C18 (ODS, Octadecyl)	1.5-10	4 x 2.0	10/pk
<b>AJO-4287</b>	C18 (ODS, Octadecyl)	1.5-10	4 x 3.0	10/pk
<b>AJO-7221</b>	C18 (ODS, Octadecyl)	1.5-10	10 x 10	3/pk
<b>AJO-7839</b>	C18 (ODS, Octadecyl)	1.5-10	15 x 21.2	ea
<b>AJO-8301</b>	C18 (ODS, Octadecyl)	1.5-10	15 x 30	ea
<b>AJO-6073</b>	C12 (Dodecyl)	1.5-10	4 x 2.0	10/pk
<b>AJO-6074</b>	C12 (Dodecyl)	1.5-10	4 x 3.0	10/pk
<b>AJO-7275</b>	C12 (Dodecyl)	1.5-10	10 x 10	3/pk
<b>AJO-7842</b>	C12 (Dodecyl)	1.5-10	15 x 21.2	ea
<b>AJO-8304</b>	C12 (Dodecyl)	1.5-10	15 x 30	ea
<b>AJO-4289</b>	C8 (Octyl, MOS)	1.5-10	4 x 2.0	10/pk
<b>AJO-4290</b>	C8 (Octyl, MOS)	1.5-10	4 x 3.0	10/pk
<b>AJO-7222</b>	C8 (Octyl, MOS)	1.5-10	10 x 10	3/pk
<b>AJO-7840</b>	C8 (Octyl, MOS)	1.5-10	15 x 21.2	ea
<b>AJO-8302</b>	C8 (Octyl, MOS)	1.5-10	15 x 30	ea
<b>AJO-4292</b>	C5 (Pentyl)	1.5-10	4 x 2.0	10/pk
<b>AJO-4293</b>	C5 (Pentyl)	1.5-10	4 x 3.0	10/pk
<b>AJO-7372</b>	C5 (Pentyl)	1.5-10	10 x 10	3/pk
<b>AJO-4298</b>	C1 (TMS)	2-9	4 x 2.0	10/pk
<b>AJO-4299</b>	C1 (TMS)	2-9	4 x 3.0	10/pk
<b>AJO-7373</b>	C1 (TMS)	2-9	10 x 10	3/pk
<b>AJO-4347</b>	Silica	—	4 x 2.0	10/pk
<b>AJO-4348</b>	Silica	—	4 x 3.0	10/pk
<b>AJO-7223</b>	Silica	—	10 x 10	3/pk
<b>AJO-7229</b>	Silica	—	15 x 21.2	ea
<b>AJO-8312</b>	Silica	—	15 x 30	ea
<b>AJO-8328</b>	HILIC	1.5-8	4 x 2.0	10/pk
<b>AJO-8329</b>	HILIC	1.5-8	4 x 3.0	10/pk
<b>AJO-8902</b>	HILIC	1.5-8	10 x 10	3/pk
<b>AJO-4301</b>	NH <sub>2</sub> (Amino, Aminopropyl)		4 x 2.0	10/pk
<b>AJO-4302</b>	NH <sub>2</sub> (Amino, Aminopropyl)		4 x 3.0	10/pk
<b>AJO-7364</b>	NH <sub>2</sub> (Amino, Aminopropyl)		10 x 10	3/pk
<b>AJO-8162</b>	NH <sub>2</sub> (Amino, Aminopropyl)		15 x 21.2	ea
<b>AJO-8309</b>	NH <sub>2</sub> (Amino, Aminopropyl)		15 x 30	ea
<b>AJO-4304</b>	CN (Cyano, Cyanopropyl)		4 x 2.0	10/pk
<b>AJO-4305</b>	CN (Cyano, Cyanopropyl)		4 x 3.0	10/pk
<b>AJO-7313</b>	CN (Cyano, Cyanopropyl)		10 x 10	3/pk
<b>AJO-8220</b>	CN (Cyano, Cyanopropyl)		15 x 21.2	ea
<b>AJO-8311</b>	CN (Cyano, Cyanopropyl)		15 x 30	ea
<b>AJO-4350</b>	Phenyl (Phenylhexyl)	1.5-10	4 x 2.0	10/pk
<b>AJO-4351</b>	Phenyl (Phenylhexyl)	1.5-10	4 x 3.0	10/pk



## Cartridges - General Purpose / Pharmaceutical

(Continued)

Part No.	Material Description	pH Stability	Dimensions L x ID(mm)	Unit
<b>AJ0-7314</b>	Phenyl (Phenylhexyl)	1.5-10	10 x 10	3/pk
<b>AJ0-7841</b>	Phenyl (Phenylhexyl)	1.5-10	15 x 21.2	ea
<b>AJ0-8303</b>	Phenyl (Phenylhexyl)	1.5-10	15 x 30	ea
<b>AJ0-8326</b>	PFP(2) (Pentafluorophenylpropyl)		4 x 2.0	10/pk
<b>AJ0-8327</b>	PFP(2) (Pentafluorophenylpropyl)		4 x 3.0	10/pk
<b>AJ0-8376</b>	PFP(2) (Pentafluorophenylpropyl)		10 x 10	3/pk
<b>AJ0-8377</b>	PFP(2) (Pentafluorophenylpropyl)		15 x 21.2	ea
<b>AJ0-8378</b>	PFP(2) (Pentafluorophenylpropyl)		15 x 30	ea
<b>AJ0-4307</b>	SCX (SA, Strong Cation Exchanger)		4 x 2.0	10/pk
<b>AJ0-4308</b>	SCX (SA, Strong Cation Exchanger)		4 x 3.0	10/pk
<b>AJ0-7369</b>	SCX (SA, Strong Cation Exchanger)		10 x 10	3/pk
<b>AJ0-8595</b>	SCX (SA, Strong Cation Exchanger)		15 x 21.2	ea
<b>AJ0-8596</b>	SCX (SA, Strong Cation Exchanger)		15 x 30	ea
<b>AJ0-4310</b>	SAX (SA, Strong Cation Exchanger)		4 x 2.0	10/pk
<b>AJ0-4311</b>	SAX (SA, Strong Cation Exchanger)		4 x 3.0	10/pk
<b>AJ0-7370</b>	SAX (SA, Strong Cation Exchanger)		10 x 10	3/pk
<b>AJ0-5808</b>	RP-1(Reversed Phase Polymer)		4 x 2.0	10/pk
<b>AJ0-5809</b>	RP-1(Reversed Phase Polymer)		4 x 3.0	10/pk
<b>AJ0-7368</b>	RP-1(Reversed Phase Polymer)		10 x 10	3/pk
<b>AJ0-8358</b>	RP-1(Reversed Phase Polymer)		15 x 21.2	ea
<b>AJ0-6075</b>	Polar-RP (Ether-linked Phenyl)		4 x 2.0	10/pk
<b>AJ0-6076</b>	Polar-RP (Ether-linked Phenyl)		4 x 3.0	10/pk
<b>AJ0-7276</b>	Polar-RP (Ether-linked Phenyl)		10 x 10	3/pk
<b>AJ0-7845</b>	Polar-RP (Ether-linked Phenyl)		15 x 21.2	ea
<b>AJ0-8307</b>	Polar-RP (Ether-linked Phenyl)		15 x 30	ea
<b>AJ0-7556</b>	Fusion-RP (C18 Polar Embedded)		4 x 2.0	10/pk
<b>AJ0-7557</b>	Fusion-RP (C18 Polar Embedded)		4 x 3.0	10/pk
<b>AJ0-7558</b>	Fusion-RP (C18 Polar Embedded)		10 x 10	3/pk
<b>AJ0-7844</b>	Fusion-RP (C18 Polar Embedded)		15 x 21.2	ea
<b>AJ0-8306</b>	Fusion-RP (C18 Polar Embedded)		15 x 30	ea
<b>AJ0-7510</b>	AQ C18 (Polar Endcapped C18)		4 x 2.0	10/pk
<b>AJ0-7511</b>	AQ C18 (Polar Endcapped C18)		4 x 3.0	10/pk
<b>AJ0-7512</b>	AQ C18 (Polar Endcapped C18)		10 x 10	3/pk
<b>AJ0-7843</b>	AQ C18 (Polar Endcapped C18)		15 x 21.2	ea
<b>AJ0-8305</b>	AQ C18 (Polar Endcapped C18)		15 x 30	ea
<b>AJ0-7596</b>	Gemini C18 (TWIN Technology)		4 x 2.0	10/pk
<b>AJ0-7597</b>	Gemini C18 (TWIN Technology)		4 x 3.0	10/pk
<b>AJ0-7598</b>	Gemini C18 (TWIN Technology)		10 x 10	3/pk
<b>AJ0-7846</b>	Gemini C18 (TWIN Technology)		15 x 21.2	ea
<b>AJ0-8308</b>	Gemini C18 (TWIN Technology)		15 x 30	ea
<b>AJ0-8367</b>	Gemini-NX (C18 TWIN-NX Technology)		4 x 2.0	10/pk
<b>AJ0-8368</b>	Gemini-NX (C18 TWIN-NX Technology)		4 x 3.0	10/pk
<b>AJ0-8369</b>	Gemini-NX (C18 TWIN-NX Technology)		10 x 10	3/pk
<b>AJ0-8370</b>	Gemini-NX (C18 TWIN-NX Technology)		15 x 21.2	ea
<b>AJ0-8371</b>	Gemini-NX (C18 TWIN-NX Technology)		15 x 30	ea
<b>AJ0-7914</b>	Gemini C6-Phenyl (TWIN Tech.)		4 x 2.0	10/pk
<b>AJ0-7915</b>	Gemini C6-Phenyl (TWIN Tech.)		4 x 3.0	10/pk
<b>AJ0-8134</b>	Oligo-RP (C18 TWIN Technology)		4 x 2.0	10/pk
<b>AJ0-8135</b>	Oligo-RP (C18 TWIN Technology)		4 x 3.0	10/pk
<b>AJ0-8136</b>	Oligo-RP (C18 TWIN Technology)		10 x 10	3/pk
<b>AJ0-8210</b>	Oligo-RP (C18 TWIN Technology)		15 x 21.2	ea
<b>AJ0-8310</b>	Oligo-RP (C18 TWIN Technology)		15 x 30	ea
<b>AJ0-8324</b>	Oligo-WAX (WA, Weak Anion Exchanger)		4 x 3.0	10/pk
<b>AJ0-8325</b>	Oligo-WAX (WA, Weak Anion Exchanger)		10 x 10	3/pk
<b>AJ0-8339</b>	Oligo-WAX (WA, Weak Anion Exchanger)		15 x 21.2	ea
<b>AJ0-8420</b>	Oligo-WAX (WA, Weak Anion Exchanger)		15 x 30	ea

(Continued on next page)

# SECURITY GUARD ORDERING INFORMATION (CONTINUED)

## Cartridges for Protein/Polypeptide Reversed Phase

For use with all silica columns for separation of proteins and peptides, such as Jupiter (Phenomenex); Vydac® 218TP, 214TP (Alltech Associates, Inc.); SynChropak® 300 C18, C4 (Eprogen, Inc.); Nucleosil® 300 Å C18, C4 (Macherey-Nagel); Hypersil® 300 Å (Thermo Hypersil-Keystone) and all other widepore or 300 Å brands.

Part No.	Material Description	pH Stability	Dimensions L x ID(mm)	Unit
<b>AJO-4320</b>	Widepore C18 (ODS)	1.5-10	4 x 2.0	10/pk
<b>AJO-4321</b>	Widepore C18 (ODS)	1.5-10	4 x 3.0	10/pk
<b>AJO-7224</b>	Widepore C18 (ODS)	1.5-10	10 x 10	3/pk
<b>AJO-7230</b>	Widepore C18 (ODS)	1.5-10	15 x 21.2	ea
<b>AJO-8313</b>	Widepore C18 (ODS)	1.5-10	15 x 30	ea
<b>AJO-4326</b>	Widepore C5 (Pentyl)	1.5-10	4 x 2.0	10/pk
<b>AJO-4327</b>	Widepore C5 (Pentyl)	1.5-10	4 x 3.0	10/pk
<b>AJO-7371</b>	Widepore C5 (Pentyl)	1.5-10	10 x 10	ea
<b>AJO-4329</b>	Widepore C4 (Butyl)	1.5-10	4 x 2.0	10/pk
<b>AJO-4330</b>	Widepore C4 (Butyl)	1.5-10	4 x 3.0	10/pk
<b>AJO-7225</b>	Widepore C4 (Butyl)	1.5-10	10 x 10	3/pk
<b>AJO-7231</b>	Widepore C4 (Butyl)	1.5-10	15 x 21.2	ea
<b>AJO-8314</b>	Widepore C4 (Butyl)	1.5-10	15 x 30	ea

## Cartridges for Silica GFC (Aqueous SEC)

For use with all silica GFC columns, such as Yarra and BioSep (Phenomenex); ZORBAX® GF-series (Agilent Technologies); Bio-Sil® (Bio-Rad)

Part No.	Material Description	pH Stability	Dimensions L x ID(mm)	Unit
<b>AJO-4487</b>	GFC-2000	2-7.5	4 x 3.0	10/pk
<b>AJO-7365</b>	GFC-2000	2-7.5	10 x 10	3/pk
<b>AJO-8588</b>	GFC-2000	2-7.5	15 x 21.2	ea
<b>AJO-4488</b>	GFC-3000	2-7.5	4 x 3.0	10/pk
<b>AJO-7366</b>	GFC-3000	2-7.5	10 x 10	3/pk
<b>AJO-8589</b>	GFC-3000	2-7.5	15 x 21.2	ea
<b>AJO-4489</b>	GFC-4000	2-7.5	4 x 3.0	10/pk
<b>AJO-7367</b>	GFC-4000	2-7.5	10 x 10	3/pk
<b>AJO-8590</b>	GFC-4000	2-7.5	15 x 21.2	ea

## Cartridges for Chiral

For use with chiral columns, such as Lux Cellulose-1, -2, -3, -4, & Amylose-2 (Phenomenex); CHIRALCEL® OD-H®, CHIRALCEL® OJ-H®, & CHIRALPAK® AD®-H (DAICEL Chemical Industries Ltd.)

Part No.	Material Description*	pH Stability	Dimensions L x ID(mm)	Unit
<b>AJO-8402</b>	Lux Cellulose-1	2-9	4 x 2.0	10/pk
<b>AJO-8403</b>	Lux Cellulose-1	2-9	4 x 3.0	10/pk
<b>AJO-8404</b>	Lux Cellulose-1	2-9	10 x 10	3/pk
<b>AJO-8405</b>	Lux Cellulose-1	2-9	15 x 21.2	ea
<b>AJO-8406</b>	Lux Cellulose-1	2-9	15 x 30	ea
<b>AJO-8398</b>	Lux Cellulose-2	2-9	4 x 2.0	10/pk
<b>AJO-8366</b>	Lux Cellulose-2	2-9	4 x 3.0	10/pk
<b>AJO-8399</b>	Lux Cellulose-2	2-9	10 x 10	3/pk
<b>AJO-8400</b>	Lux Cellulose-2	2-9	15 x 21.2	ea
<b>AJO-8401</b>	Lux Cellulose-2	2-9	15 x 30	ea

(Continued on next page)

## Cartridges for Chiral (cont'd)

Part No.	Material Description*	pH Stability	Dimensions L x ID(mm)	Unit
AJ0-8621	Lux Cellulose-3	2-9	4 x 2.0	10/pk
AJ0-8622	Lux Cellulose-3	2-9	4 x 3.0	10/pk
AJ0-8623	Lux Cellulose-3	2-9	10 x 10.0	3/pk
AJ0-8624	Lux Cellulose-3	2-9	15 x 21.2	ea
AJ0-8625	Lux Cellulose-3	2-9	15 x 30.0	ea
AJ0-8626	Lux Cellulose-4	2-9	4 x 2.0	10/pk
AJ0-8627	Lux Cellulose-4	2-9	4 x 3.0	10/pk
AJ0-8628	Lux Cellulose-4	2-9	10 x 10.0	3/pk
AJ0-8629	Lux Cellulose-4	2-9	15 x 21.2	ea
AJ0-8630	Lux Cellulose-4	2-9	15 x 30.0	ea
AJ0-8471	Lux Amylose-2	2-9	4 x 2.0	10/pk
AJ0-8470	Lux Amylose-2	2-9	4 x 3.0	10/pk
AJ0-8472	Lux Amylose-2	2-9	10 x 10	3/pk
AJ0-8473	Lux Amylose-2	2-9	15 x 21.2	ea
AJ0-8474	Lux Amylose-2	2-9	15 x 30	ea

\* Lux Cellulose-1 is cellulose tris(3,5-dimethylphenylcarbamate)  
Lux Cellulose-2 is cellulose tris(3-chloro-4-methylphenylcarbamate)  
Lux Cellulose-3 is cellulose tris(4-methylbenzoate)  
Lux Cellulose-4 is cellulose tris(4-chloro-3-methylphenylcarbamate)  
Lux Amylose-2 is amylose tris(5-chloro-2-methylphenylcarbamate)

## Cartridges for Carbohydrate / Organic Acid

For Organic acid and carbohydrate analysis, such as Rezex™ (Phenomenex); Aminex® (Bio-Rad); Sugar-Pak™ (Waters).

Part No.	Material Description	pH Stability	Dimensions L x ID(mm)	Unit
AJ0-4490	Carbo-H <sup>+</sup>	1 - 8	4 x 3.0	10/pk
AJ0-8888	Carbo H <sup>+</sup>	1 - 8	15 x 21.2	ea
AJ0-4491	Carbo-Ag <sup>+</sup> *	Neutral	4 x 3.0	10/pk
AJ0-8891	Carbo Ag <sup>+</sup> *	Neutral	15 x 21.2	ea
AJ0-4492	Carbo-Pb <sup>+2</sup>	Neutral	4 x 3.0	10/pk
AJ0-8890	Carbo Pb <sup>+2</sup>	Neutral	15 x 21.2	ea
AJ0-4493	Carbo-Ca <sup>+2</sup>	Neutral	4 x 3.0	10/pk
AJ0-8889	Carbo Ca <sup>+2</sup>	Neutral	15 x 21.2	ea

\*For use with saccharide and oligosaccharide columns in Ag<sup>+</sup> form.

## Replacement Parts

Part No.	Description	Unit
AJ0-4283	PEEK Ferrules	3/pk
AJ0-4285	Stacking Rings	2/pk
AQ0-1389	PEEK Fingertight Male Nuts	10/pk
AJ0-4284	Security Guard Wrenches	2/pk
AQ0-8374	PREP Coupler, SS w/ PEEK Ferrule Inserts 10-32 Threads, 1/16 in. OD x 0.020 in. ID	ea
AQ0-8375	Replacement Ferrule Inserts, for PREP Coupler, PEEK, 0.020 in. ID	10/pk
AQ0-8222	PREP Replacement O-Rings, Kalrez® For 15 x 21.2 mm SG Holder, Size 2-021	2/pk
AQ0-8318	PREP Replacement O-Rings, Kalrez® For 15 x 30 mm SG Holder, Size 2-025	2/pk

## High Pressure Applications

For applications where pressures exceed 3500 psi (241bar), or for use with core-shell, non-porous, or < 3 µm fully porous media, please choose SecurityGuard ULTRA. Contact Phenomenex or your local Phenomenex representative for more information.


## HPLC SYSTEM TEST KIT



- Diagnose hardware problems rapidly and easily
- Avoid unnecessary and costly system repairs
- Convenient benchmark testing of HPLC systems using a C18 column standard
- Test system setup and hardware connections
- Quickly isolate method development problems
- Reduce instrument downtime

Each kit contains the following:

1. Phenomenex 5  $\mu$ m C18, 50 x 4.6 mm HPLC column
2. Five vials of Isocratic Test Mix
3. Five vials of Gradient Test Mix



## ORDERING INFORMATION

Part No.	Description	Unit
CHO-1684	HPLC System Test Kit, Reversed Phase, includes: C18 column, isocratic and gradient test mixes	ea
CHO-1685	Isocratic Test Mix	5/pk
CHO-1686	Gradient Test Mix	5/pk

## COLUMN PERFORMANCE CHECK STANDARDS



- Convenient way to check column performance
- Affordable and easy to use

Phenomenex offers a comprehensive line of column performance check standards to help you evaluate column performance. We recommend using the check standards to verify performance of all columns upon receiving them and periodically over the lifetime of the column. Test conditions are located in the column jacket.

### NORMAL PHASE

Part No. **ALO-3033**

(For Si, NH<sub>2</sub>, NO<sub>2</sub>, Alumina, PAC, and Luna CN)

**Unit quantity:** 2 mL

**Contains:** Meta-xylene, Nitrobenzene

### REVERSED PHASE 1

Part No. **ALO-3034**

(For C1, C18, CN and Phenyl)

**Unit quantity:** 2 mL

**Contains:** Uracil, Benzamide, Benzophenone, Biphenyl

### REVERSED PHASE 2

Part No. **ALO-3045**

(For Prodigy C8, ODS(2), ODS(3); Luna C5, C8, C18, Phenyl-Hexyl, PFP(2); Jupiter C4, C5, C18, Proteo; Columbus C8, C18; Aqua; Synergi; PhenoSphere-NEXT C8, C18; Gemini C18, C6-Phenyl; Gemini NX-C18; Clarity Oligo-RP, Oligo-MS; Kinetex C8, XB-C18, Phenyl-Hexyl, C18, PFP; Aeris PEPTIDE XB-C18; 4.6 mm ID Aeris WIDEBORE XB-C18, XB-C8, C4)

**Unit quantity:** 2 mL

**Contains:** Uracil, Acetophenone, Toluene, Naphthalene  
(Please refer to the QC Test Data for specific test conditions for Jupiter, Aeris, Kinetex, and Luna)

---

## COLUMN PERFORMANCE CHECK STANDARDS (CONTINUED)

---

### **AERIS NARROW ID** Part No. **ALO-8931**

(For 2.1 mm ID Aeris WIDEPORE XB-C18, XB-C8, C4)

---

**Unit quantity:** 2 mL

**Contains:** Uracil; Acetophenone; Toluene; Naphthalene;  
Acenaphthalene  
(2.5 mg/mL)

---

### **HILIC PHASE** Part No. **ALO-8317**

(For Luna HILIC; Kinetex HILIC)

---

**Unit quantity:** 2 mL

**Contains:** Toluene, Uracil, Cytosine

---

### **CARBOHYDRATE MIX 1** Part No. **ALO-3035**

(For Rezex RNM, RAM and other carbohydrate analysis columns)

---

**Unit quantity:** 2 mL

**Contains:** Maltotriose Hydrate, Maltose, Ribitol

---

### **CARBOHYDRATE MIX 2** Part No. **ALO-3036**

(For Rezex RPM and other carbohydrate analysis columns)

---

**Unit quantity:** 2 mL

**Contains:** Melezitose, Glucose, Fructose, Ribitol

---

### **CARBOHYDRATE MIX 3** Part No. **ALO-3037**

(For Rezex RCM, RCU, and other carbohydrate analysis columns)

---

**Unit quantity:** 2 mL

**Contains:** Melezitose, Maltose, Glucose, Mannose,  
Fructose, Ribitol

---

### **OLIGOSACCHARIDE STANDARD** Part No. **ALO-3038**

(For Rezex RSO, RNO, and other oligosaccharide analysis columns)

---

**Unit quantity:** 2 mL

**Contains:** Light corn syrup

---

### **ORGANIC ACID STANDARD** Part No. **ALO-3039**

(For Rezex ROA and other organic acid analysis columns)

---

**Unit quantity:** 2 mL

**Contains:** Oxalic acid, Succinic acid, Citric acid,  
Formic acid, Tartaric acid, Acetic acid

---

### **CATION-EXCHANGE** Part No. **ALO-3040**

(For SCX, SA, CM)

---

**Unit quantity:** 2 mL

**Contains:** Uracil, Cytosine

---

### **ANION-EXCHANGE** Part No. **ALO-3041**

(For SAX, SB, DEAE, PEI)

---

**Unit quantity:** 2 mL

**Contains:** Uridine, UMP

---

# COLUMN PERFORMANCE CHECK STANDARDS (CONTINUED)

---

## AQUEOUS SEC 1

Part No. **ALO-3042**

(For Yarra SEC, BioSep-SEC-S, and other protein SEC columns)

---

**Unit quantity:** Dry; Reconstituted to 2 mL  
**Contains:** Bovine thyroglobulin  
Human gamma globulin (contains IgA and IgG)  
Ovalbumin  
Myoglobin  
Uridine

*(reconstitute with 1 mL of 100 mM Sodium Phosphate pH 6.8)*

---

## AQUEOUS SEC 2

Part No. **ALO-3043**

(For PolySep GFC-P and other aqueous-soluble analysis columns)

---

**Unit quantity:** 2 mL  
**Contains:** Ethylene Glycol

---

## STAR-ION A300

Part No. **ALO-3420**

---

**Unit quantity:** 2 mL  
**Contains:**

<u>Conc. (mg/mL)</u>			
	Nitrite		20
Fluoride	5	Sulfate	20
Nitrate	20	Bromide	20
Chloride	10		
Phosphate	30		

---

## POLYMERX RP-1

Part No. **ALO-7260**

---

**Unit quantity:** 2 mL  
**Contains:**

<u>Conc. (mg/mL)</u>	
Cytosine	13
Uracil	13
Uridine	33

---

## ONYX MONOLITHIC REVERSED PHASE

Part No. **ALO-7836**

---

**Unit quantity:** 2 mL  
**Contains:**

<u>Conc. (µg/mL)</u>	
Thiourea	10
Progesterone	100
Anthracene	10

---

## ONYX MONOLITHIC NORMAL PHASE

Part No. **ALO-7835**

---

**Unit quantity:** 2 mL  
**Contains:**

<u>Conc. (µg/mL)</u>	
Toluene	21.75
Nitrobenzene	150.00
2-Nitroanisol	0.18

---

## COLUMN PERFORMANCE CHECK STANDARDS (CONTINUED)

### CHIRAL TEST MIX 1

Part No. **ALO-3046**

Applicable to the following Chirex columns:  
3001, 3005

**Unit quantity:** 2 mL

**Contains:**

1. S-(+)-2,2,2-trifluoro-1-(9-anthryl)  
ethanol CAS [60646-30-2]
2. R-(-)-2,2,2-trifluoro-1-(9-anthryl)  
ethanol CAS [53531-34-3]

### CHIRAL TEST MIX 2

Part No. **ALO-3047**

Applicable to the following Chirex columns:  
3010, 3011, 3012

**Unit quantity:** 2 mL

**Contains:** N-dansyl-DL-valine  
(cyclohexylammonium salt)  
CAS[84540-67-0]

### CHIRAL TEST MIX 3

Part No. **ALO-3048**

Applicable to the following Chirex columns:  
3014, 3017, 3018, 3019, 3020, 3022

**Unit quantity:** 2 mL

**Contains:**

1. (R)-(-)-N-(3,5-Dinitrobenzoyl)- $\alpha$ -  
methylbenzylamine CAS [69632-32-2]
2. (S)-(-)-N-(3,5-Dinitrobenzoyl)- $\alpha$ -  
methylbenzylamine CAS[69632-31-1]

### CHIRAL TEST MIX 4

Part No. **ALO-3049**

Applicable to the following Chirex column:  
3126

**Unit quantity:** 2 mL

**Contains:** DL-Aspartic Acid CAS [617-45-8]

### CHIRAL TEST MIX 5

Part No. **ALO-8412**

Applicable to the following Lux columns:  
Lux Cellulose -1,-2,-3,-4, Lux Amylose-2

**Unit quantity:** 2 mL

**Contains:** trans-Stilbene oxide CAS [1439-07-2]

## PART XIII - VIALS

### VEREX VIAL PRODUCTS CERTIFIED VIALS, CAPS, SEPTA, AND INSERTS

#### CERTIFIED FOR DEMANDING METHODOLOGIES:

- Regulated Methods
- High Sensitivity LC and GC
- Mass Spectrometry




#### ORDERING INFORMATION

#### 12 x 32 mm, 11 mm Crimp-Top Vials and Closures

##### CRIMP-TOP VIALS, 2.0 mL

- Cleaner vials eliminate ghost peaks and contaminants
- Used with most autosamplers, including Agilent®, Thermo Scientific®, Waters®
- Larger-opening “wide-mouth” style prevents broken needles and system downtime
- Precision neck improves crimping



Description	1000/pk
<b>Standard Opening</b>	
Vial, Crimp, 2 mL Clear, No Patch	AR0-3700-13
Vial, Crimp, 2 mL Clear, w/ Patch	AR0-3710-13
Vial, Crimp, 2 mL Amber, w/ Patch	AR0-3711-13
<b>Wide Mouth Opening</b>	
Vial, Crimp, 2 mL Wide Mouth, Clear, No Patch	AR0-37K0-13
Vial, Crimp, 2 mL Wide Mouth, Clear, w/ Patch	AR0-37L0-13
Vial, Crimp, 2 mL Wide Mouth, Amber, No Patch	AR0-37K1-13
Vial, Crimp, 2 mL Wide Mouth, Amber, w/ Patch	AR0-37L1-13

##### SEALS / CLOSURES FOR CRIMP-TOP VIALS

- Excellent for volatile samples
- Extra clean to eliminate contamination
- Colored aluminum



Description	1000/pk
Seal, 11mm Diameter, Crimp, PTFE/Silicone, silver	AR0-5780-13
Seal, 11mm Diameter, Crimp, PTFE/Silicone/PTFE, silver	AR0-5760-13
Seal, 11mm Diameter, Crimp, PTFE/Rubber, silver	AR0-5740-13
Seal, 11mm Diameter, Crimp, PTFE/Rubber, blue	AR0-5742-13
Seal, 11mm Diameter, Crimp, PTFE/Rubber, red	AR0-5741-13
Seal, 11mm Diameter, Crimp, PTFE/Rubber, green	AR0-5743-13
Seal, 11mm Diameter, Crimp, PTFE/Rubber, gold	AR0-574G-13
Seal, 11mm Diameter, Crimp, PTFE, silver	AR0-5710-13



Need help matching your current vials and caps to Verex? Visit [www.phenomenex.com/VialMatch](http://www.phenomenex.com/VialMatch)



## 12 x 32 mm, 9-425 (9 mm) Screw-Top Vials and Caps

### 9-425 SCREW-TOP VIALS, 2.0 mL

- Used with most autosamplers, including Agilent®, Thermo Scientific®, Waters® and many others
- Performs as well as crimp or snap vials
- Offers improved cap convenience and accessibility (easy on, easy off)



Description	1000/pk
Vial, 9 mm Screw, 2 mL Clear, No Patch	ARO-3900-13
Vial, 9 mm Screw, 2 mL Amber, No Patch	ARO-3901-13
Vial, 9 mm Screw, 2 mL Clear, w/ Patch	ARO-3910-13
Vial, 9 mm Screw, 2 mL Amber, w/ Patch	ARO-3911-13
Vial, 9 mm Screw, 2 mL Clear, w/ Patch, Silanized	ARO-3960-13

### BONDED-IN CAPS FOR 9-425 SCREW-TOP VIALS

- Bonded septa caps eliminate costly liner/septa fallout
- Prevents rework and wasted productivity with perfect-fit septa
- Saves instrument downtime



Description	1000/pk
Cap (pre-assembled), 9 mm, w/ Bonded-in PTFE/Silicone septa, black	ARO-8957-13-B
Cap (pre-assembled), 9 mm, w/ Bonded-in PTFE/Silicone septa, blue	ARO-8952-13-B
Cap (pre-assembled), 9 mm, w/ Bonded-in PTFE/Silicone septa, natural	ARO-8956-13-B
Cap (pre-assembled), 9 mm, w/ Bonded-in PTFE/Silicone septa, red	ARO-8951-13-B
Cap (pre-assembled), 9 mm, w/ Bonded-in PTFE/Silicone preSlit septa, black	ARO-8977-13-B
Cap (pre-assembled), 9 mm, w/ Bonded-in PTFE/Silicone preSlit septa, blue	ARO-8972-13-B

## ADDITIONAL CERTIFIED VIALS AND CAPS

### ALSO AVAILABLE:

- 12 x 32 mm, 2 mL Vials
  - 11 mm limited volume crimp-top vials
  - 11 mm limited volume screw-top vials
  - 11 mm snap-top vials
  - 11 mm limited volume snap-top vials
  - 8 mm screw-top vials
  - 10 mm screw-top vials
  - 10 mm limited volume screw-top vials
  - 13 mm screw-top vials
  - 13 mm limited volume screw-top vials
- VOA / ASE assembled vial kits and storage vial
- Headspace vials
- Plastic vials
- Shell vials
- Vial inserts



**Learn More.** Visit [Phenomenex.com/Verex](http://Phenomenex.com/Verex) for additional product selection and detailed information.

## PART XIV - SOLID PHASE EXTRACTION (SPE)

Increase column and instrument life by injecting samples cleaned-up with Strata®.

### STRATA™-X Polymeric Sorbents

Tubes and 96-Well Plates

- Deconditioning Resistant
- Low Elution Volumes
- High Analyte Capacity



Strata™-X and -XL for simplified cleanup of polar and non-polar compounds

Strata™-X-C and -XL-C for selective extraction of basic compounds

Strata™-X-CW and -XL-CW for bases (including quaternary amines)

Strata™-X-A and -XL-A for cleanup of weak acids

Strata™-X-AW -XL-AW for acids

Strata™-X-Drug B for basic drugs of abuse

Strata™-X-Drug N for neutral drugs of abuse

### STRATA® Traditional Sorbents

Tubes and 96-Well Plates

- Optimal Flow
- Lot-to-Lot Reproducibility
- Wide Range of Selectivity
- Available chemistries include:  
C18-E, C18-U, C18-T, C8,  
Phenyl, SDB-L, CN, Si-1, WCX, FI-PR, NH<sub>2</sub>, SAX, SCX, Melamine



### STRATA® Flash Sorbents

- Polar & Non-polar Phases
- Narrow Particle Range Distribution
- Can be used for Direct Scale-up



Strata® Giga™ Tubes available in

12, 20, 60 & 150 mL formats

Sepra™ Bulk available in gram to multi-kilogram quantities

### STRATA® On-line Cartridges



- Rapid Extraction and Concentration
- Direct Inject Analysis
- Easily Automated



Strata™-X for polar and non-polar compounds

Strata™-X-C for weak bases

Strata™-X-CW for strong bases

Strata® C18 for non-polar compounds

Strata® C8 for compounds of intermediate polarity

**Subject to Phenomenex Standard Terms & Conditions which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions)** Gemini, Jupiter, Luna, Lux, Kinetex and Strata are registered trademarks of Phenomenex. Aqua, Axia, BioSep, Columbus, Giga, Aeris, Onyx, Phenex, Phenogel, PhenoSphere, PolymerX, PolySep, Prodigy, Rezex, SecurityGuard, Sepra, Star-Ion, Strata-X, Synergi, ThermoSphere, Verex, Yarra, TWIN-NX and TWIN Technology are trademarks of Phenomenex. Alltima is a trademark of Alltech Associates, Inc. Bio-Sil and Aminex are registered trademarks of Bio-Rad Laboratories. Chirex is a trademark of Chirex, Inc. licensed to Phenomenex. Cyrolite is a registered trademark of CY/RO Industries. Hypersil is a registered trademark of Thermo Hypersil-Keystone. Nucleosil is a registered trademark of Macherey-Nagel. Sugar-Pak is a trademark of Waters Corporation. Symmetry is a registered trademark of Waters Corporation. SynChropak is a registered trademark of Eprogen, Inc. Teflon is a registered trademark of E.I. du Pont de Nemours and Co. Vydac is a registered trademark of Alltech Associates, Inc. ZORBAX is a registered trademark of Agilent Technologies. CHIRALCEL, OD-H, OJ-H, CHIRALPAK, and AD are registered trademarks of DAICEL Chemical Industries, Ltd. Onyx is a product based on monolithic technology under license from Merck KGaA, Darmstadt, Germany.

© 2013 Phenomenex. All rights reserved.

## PART XV - HPLC ACCESSORIES

### ACCESSORIES

- Backpressure Regulators
- Biocompatible / Metal-free products
- Connectors and Splitters
- Filtration Products
- Injectors and Injector Loops
- Membrane Filters
- Mobile Phase Handling Devices
- Polymer Calibration Standards / Kits
- Rotor Seals, Stators, etc.
- Solvent Reservoir and Reagent Bottles
- SPE Consumables, Tube & Plate Manifolds
- Switching Valves
- Syringes
- Syringe Filters
- Tools
- Tubing, Fittings, Frits and Unions
- Valves (Injection, Switching)
- Vials, Caps and Septa

---


### EQUIPMENT

- Column Chiller-Heater
- Column Heater
- Degasser

*For ordering and additional information, please contact your Phenomenex Technical Consultant.*

#### SINGLE COLUMN HEATER THERMASPHERE™ TS-130

- Compact, low-cost heater precisely controls temperature from 25-90 °C
- Improves reproducibility and chromatographic results
- Reduces analyte identification errors
- Improves baseline and overall detector performance
- Improves peak efficiency and analyte quantitation (especially at low levels)
- Improves the ruggedness of the separation (within-lab and lab-to-lab)



---

### ORDERING INFORMATION

#### ThermaSphere™ TS-130

Part No.	Description
<b>EH0-7057</b>	ThermaSphere TS-130 HPLC Column Heater 25-90 °C, 95 to 265 VAC, 50/60 Hz
<b>EH0-7058</b>	Stand for ThermaSphere TS-130 HPLC Column Heater

*More Accessories available. See Phenomenex Catalog for details.*

# PHENOMENEX WARRANTY

Phenomenex products are warranted to meet the stated performance and quality and to be free of defects in material and workmanship. They are not warranted, nor does Phenomenex assume liability, if misused. **NO OTHER WARRANTY OR REPRESENTATION IS IMPLIED OR EXPRESSED BY PHENOMENEX FOR ITS PRODUCTS WITH RESPECT TO MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR ANY OTHER MATTER. PHENOMENEX SHALL NOT UNDER ANY CIRCUMSTANCES BE LIABLE FOR ANY INCIDENTAL, CONSEQUENTIAL, OR COMPENSATORY DAMAGE ARISING FROM THE USE OF, OR IN CONJUNCTION WITH, ITS PRODUCTS.**

The maximum liability which can be assumed by Phenomenex for breach of warranty shall be the invoice price of the product.

---

## SPECIFIC WARRANTIES ON HPLC COLUMNS

Phenomenex warrants its quality columns in accordance with the following terms and conditions. Phenomenex will repack, replace, or refund charges on any column (at our discretion), at no cost if a column fails to perform satisfactorily. Columns being returned must have prior return authorization granted by Phenomenex. Defective products must be accompanied by a written explanation of failure. Approval is subject to the following exclusions:

- All columns must be tested upon receipt and all deficiencies must be reported to Phenomenex no later than 15 days after the date of receipt of the column.
- Maximum warranty period is limited to 90 days on HPLC columns unless previously agreed upon. However, **COLUMNS MAY NOT BE RETURNED FOR REFUND OR CREDIT AFTER 45 DAYS AND WITHOUT PRIOR AUTHORIZATION.**
- Removal of column end-fittings automatically voids column warranty.
- Column performance warranty is limited to the conditions of the original test chromatograms.
- Physical damage to the column due to misuse, abuse, or mishap, including mechanical shock.
- Chemical damage to the packing material due to operation at incorrect chemical conditions, temperatures, or pressures.
- Failure due to high backpressures caused by improper solvent or sample filtration practices causing particulate build-up or precipitation in the column or end-fitting.
- Incorrect selection of packing material made by customer for their particular use or incompatibility of equipment, etc.
- For products supplied by but not manufactured by Phenomenex, the warranty is limited by the terms of the original manufacturer's warranty.





**www.phenomenex.com**

Phenomenex products are available worldwide. For the distributor in your country,  
contact Phenomenex USA, International Department by telephone, fax or email: [international@phenomenex.com](mailto:international@phenomenex.com).

QUALITY MANAGEMENT SYSTEM

CERTIFIED BY DNV

ISO 9001:2008

 **phenomenex**<sup>®</sup>  
*...breaking with tradition<sup>SM</sup>*

[www.phenomenex.com](http://www.phenomenex.com)

[support@phenomenex.com](mailto:support@phenomenex.com)