KINETEX®

**Column Care and
Optimization Notes**

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SPECIFICATIONS

Shipping Solvent

- Kinetex C18, EVO C18, XB-C18, C8, Biphenyl, Phenyl-Hexyl, F5 and PFP columns are shipped in Acetonitrile / Water ($\geq 50:50$ v/v; exact composition dependent on column dimensions)
- Kinetex HILIC is shipped in Acetonitrile / Water 90:10 (v/v)

Test Certificate

Each Kinetex C18, EVO C18, XB-C18, C8, Biphenyl, Phenyl-Hexyl, F5 and PFP column is individually tested before shipment. A test certificate showing the separation parameters for the reversed phase test mixture containing uracil, acetophenone, toluene, and naphthalene can be found online at www.phenomenex.com/mysupport.

Each Kinetex HILIC column is individually tested before shipment. A test certificate showing the separation parameters for the HILIC test mixture containing toluene, uracil, and cytosine can be found online at www.phenomenex.com/mysupport.

If the performance of your Kinetex column is not similar to the test certificate, please review the system optimization tips in this guide or contact Phenomenex.

OPERATING PARAMETERS

Column Installation

The arrows on the column tag indicate the flow direction. Do not operate or back-flush a Kinetex 1.3 μm column in the opposite direction to the one indicated by the arrow on the column tag. Phenomenex recommends the use of HPLC / UHPLC Sure-Lok™ High Pressure PEEK male nut fittings for installation of Kinetex columns on HPLC / UHPLC systems. The convenient one-piece design (AQO-8503) is pressure rated to 12,000 psi (827 bar), while the 3-piece unit (AQO-8504 with AQO-8505), which contains a PEEK nut, a ferrule, and stainless steel gripping ring, will provide leak-free connections up to 19,000 psi (1,310 bar). A handy fitting tightening tool (AQO-8530) is available to facilitate achievement of a leak-free connection.

Mobile Phase Restrictions

- Kinetex EVO C18 columns are stable from pH 1 - 12 and can be used with typical reversed phase mobile phase (aqueous methanol, aqueous acetonitrile or appropriate aqueous buffer / methanol or aqueous buffer / acetonitrile mixtures)
- Kinetex C18, XB-C18, C8, Biphenyl and Phenyl-Hexyl columns are stable from pH 1.5 to 10* and can be used with typical reversed phase mobile phase (aqueous methanol, aqueous acetonitrile or appropriate aqueous buffer / methanol or aqueous buffer / acetonitrile mixtures)
- Kinetex F5 and PFP columns are stable from pH 1.5 to 8.5 and can be used with typical reversed phase mobile phase (aqueous methanol, aqueous acetonitrile or appropriate aqueous buffer / methanol or aqueous buffer / acetonitrile mixtures)
- Kinetex HILIC columns are stable from pH 2.0-7.5 and can be used with typical HILIC mobile phase (acetonitrile / aqueous buffer mixtures)

* For isocratic conditions only.

Use only high purity reagents and high quality chromatography grade solvents to prepare mobile phase. Trace impurities can dramatically degrade column lifetime. Degas and filter all mobile phase prior to use. Ensure sample (and matrix) are completely soluble / miscible with mobile phase. Immiscible solvents or buffer salt precipitation can permanently damage the column.

Kinetex® Core-Shell Technology columns provide performance gains on any LC system.

Please carefully read these care and use notes to obtain the best results on your system.

Avoid:

- Operating below the lower pH limit specified for each Kinetex column as this will hydrolyze the bonded phase
- Operating above the upper pH limit specified for each Kinetex column as this will dissolve the silica
- Operating a Kinetex 1.3 μm column in the opposite direction to what is indicated by the arrow on the column tag
- Immiscible solvents and buffers
- Sudden pressure changes

Operating Backpressure

The maximum operating pressure for Kinetex 5 μm and 2.6 μm columns is 8,700 psi (600 bar)*, and for Kinetex 1.7 μm and 1.3 μm columns it is 15,000 psi (1,000 bar). The mobile phase flow rate should be set such that the column backpressure does not exceed this maximum operating pressure. Note that operating at or near the maximum pressure will result in shorter column lifetimes.

* 2.1 mm ID columns are pressure stable up to 1,000 bar.

Operating Temperatures

Kinetex columns can be used up to a maximum temperature of 60 °C. Reduced mobile phase viscosity (and backpressure) and increased mass transfer rates will be realized above 25 °C. When operating at high pH (> 8), lower operating temperatures are recommended for longer column lifetime. Note that operating at or near the maximum temperature will result in shorter column lifetimes.

Column Cleaning and Regeneration

If an increase in operating backpressure is observed, reverse flush the column (do not try this on any Kinetex 1.3 μm column or any other manufacturers' columns) with reduced flow rates indicated below:

Column ID	Flow Rate
2.1 mm	0.1 mL/min
3.0 mm	0.3 mL/min
4.6 mm	0.5 mL/min

Reversed Phase (C18, EVO C18, XB-C18, C8, Biphenyl, Phenyl-Hexyl, F5 and PFP)

Kinetex reversed phase columns can be cleaned by flushing with 10 to 20 column volumes of the following solvent mixtures in succession:

- 1) 5:95 Acetonitrile / Water (or Methanol / Water) for buffer removal
- 2) 95:5 Acetonitrile / Water (or Methanol / Water)
- 3) THF (Tetrahydrofuran)
- 4) 95:5 Acetonitrile / Water (or Methanol / Water)
- 5) 5:95 Acetonitrile / Water (or Methanol / Water)
- 6) Equilibrate with mobile phase

HILIC

Kinetex HILIC columns can be cleaned by flushing with 10 to 20 column volumes of:

- 1) 95 % Water / 5 % Acetonitrile (For buffer removal)
- 2) 95 % 100 mM Ammonium Acetate, pH 5.8 / 5 % Acetonitrile
- 3) 95 % Water / 5 % Acetonitrile
- 4) Equilibrate with mobile phase

Column Storage

Reversed Phase (C18, EVO C18, XB-C18, C8, Biphenyl, Phenyl-Hexyl, F5 and PFP)

Column storage for a period longer than several days is recommended in $\geq 50\%$ (v/v) acetonitrile or methanol in water. If the mobile phase contained buffer salt, then flush the column with 10 to 20 column volumes of water / acetonitrile or water / methanol to remove the buffer salts before storage. After flushing the column, insert column end plugs securely to prevent evaporation and drying out of the column bed.

HILIC

Column storage for a period longer than several days is recommended in 90% (v/v) acetonitrile in water. If the mobile phase contained buffer salt, then flush the column with 10 to 20 column volumes of 90:10 acetonitrile / water to remove the buffer salts before storage. After flushing the column, insert column end plugs securely to prevent evaporation and drying out of the column bed.

Kinetex 1.3 μm Operating Tips

Every care should be made to prevent and eliminate microbial growth in mobile phase reservoirs and UHPLC system plumbing when using a Kinetex 1.3 μm column. Some general recommended tips include:

Solvents: Use ultra-pure, high quality UHPLC solvents

Aqueous Mobile Phases: Change daily or add 5-10% organic modifier to water

Mobile Phase-Reservoirs: Rinse aqueous reservoirs with methanol prior to refilling

UHPLC System: Flush with isopropanol or methanol weekly

OPTIMIZING THE PERFORMANCE OF YOUR KINETEX® COLUMN

Flow Rate

Kinetex columns are capable of maintaining high efficiencies with increasing flow rate. To reduce analysis times using Kinetex use shorter length columns (30, 50 or 75 mm) and increase the mobile phase flow rate.

There is an optimal Kinetex column for your system and operating conditions.

Visit

www.phenomenex.com/optimize

to determine the starting Kinetex column dimensions for your method.

System Optimization

Kinetex® 2.6 µm core-shell columns operate comfortably within the pressure limits of conventional HPLC instruments and rival the performance of sub-2 µm fully porous particles achieved on UHPLC instruments. However, to maximize the benefits from your Kinetex core-shell column do the following:

- Minimize sample dispersion before the column
 1. Gradient elution will minimize dispersion
 2. In isocratic elution, use an injection solvent that's weaker than your mobile phase
- Optimize detector settings by adjusting the scan rate and/or the time constant to the fastest possible setting such that signal-to-noise (s/n) is not adversely affected.
- Minimize the extra-column volume between the injector and the column, and between the column and the detector
 1. Minimize the length of all connection tubing
 2. Use 0.12 mm ID (0.005 in.) tubing whenever possible
 - a) 0.17 mm ID (0.007 in.) tubing is acceptable
 - b) Avoid, if possible, using 0.25 mm ID (0.010 in.) tubing
 3. Use extremely low dead-volume fittings
 4. Ensure all connecting tubing is seated properly at every connection
- Use a micro volume flow-cell
 1. Standard flow cells on conventional LC systems can be >10 µL
 2. For best results, replace standard flow cells with <3 µL flow cells (<2 µL when using 2.1 mm ID columns)

NOTE: To convert existing HPLC and UHPLC methods for use on Kinetex columns we have created a conversion calculator, available on the Phenomenex website: www.phenomenex.com/optimize

Extending Column Lifetime

Phenomenex recommends the use of the SecurityGuard™ ULTRA guard cartridge system to extend the lifetime of your Kinetex column, especially with samples extracted from complex matrixes. Ideally, samples must be completely dissolved in the mobile phase or filtered through a syringe filter of approximately 0.20 µm porosity.

SecurityGuard ULTRA Cartridge Holder

Ordering Information

Part No.	Description	Unit
AJ0-9000	SecurityGuard ULTRA Cartridge Holder	ea

Cartridge Holder

SecurityGuard ULTRA Cartridges

Ordering Information

Cartridges (3pk)	Column ID (mm)		
	2.1	3.0	4.6
EVO C18	AJ0-9298	AJ0-9297	AJ0-9296
C18	AJ0-8782	AJ0-8775	AJ0-8768
C8	AJ0-8784	AJ0-8777	AJ0-8770
PFP	AJ0-8787	AJ0-8780	AJ0-8773
F5	AJ0-9322	AJ0-9321	AJ0-9320
Biphenyl	AJ0-9209	AJ0-9208	AJ0-9207
HILIC	AJ0-8786	AJ0-8779	AJ0-8772
Phenyl	AJ0-8788	AJ0-8781	AJ0-8774

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions which may be viewed at www.Phenomenex.com/TermsAndConditions

Trademarks

Kinetex is a registered trademark, SecurityGuard is a trademark of Phenomenex. Kinetex EVO is patented by Phenomenex. U.S. Patent Nos. 7,563,367 and 8,658,038 and foreign counterparts

 Phenomenex products are available worldwide.

www.phenomenex.com/mysupport

© 2015 Phenomenex, Inc. All rights reserved.

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www.phenomenex.com/mysupport

TAKE ADVANTAGE OF OUR KINETEX® WEB SERVICES TODAY!

- Product information
- Applications and Technical notes
- Installation and Optimization videos
- Productivity and Method Optimization calculators
- Quote request
- And so much more...

The screenshot displays the Phenomenex website interface. At the top, the Phenomenex logo is on the left, and navigation links for 'My Dashboard' and 'Quick Order' are on the right. A search bar is positioned below the navigation. The main header features the text 'Shockingly Better Performance Than Your Current LC Column!' and the Kinetex logo. Below this, a horizontal menu includes 'Easy Column Selection', 'Core-Shell Technology', 'Technical Resources', 'Application Search', 'Ask Ball-Guy', 'Care Club', and 'Order Now'. The central content area is titled 'Column Installation and Connectivity' and contains four video thumbnails with captions: 'Kinetex UHPLC Column Installation into a Waters AGILITY II', 'Kinetex Column Installation with Mist Tubing & Inlet Fittings', 'Kinetex Column Installation with PEEK Tubing & PEEK FINGERPOINT™ Fittings', and 'Kinetex HPLC Column Installation into a Shimadzu B Phenomenex B-LC-30A'. A sidebar on the left lists various technical resources. At the bottom, there are three promotional boxes for 'What's New!', 'Kinetex columns', and 'Product Videos'. The footer contains a site map, contact information, and social media icons.

Want to witness the power of Kinetex columns at your company?

Request an on-site demonstration or seminar at:

www.phenomenex.com/KinetexDeminar

The screenshot shows the Phenomenex website interface. At the top left is the Phenomenex logo with the tagline "...making with nature". To the right are links for "Log In", "My Account", "Quick Order", and a shopping cart icon. Below the logo is a navigation menu with links for "Products", "Industries", "News and Events", "Application & Support", and "Our Company". A search bar is located in the top right. The main banner features the Kinetex logo and the headline "Shockingly Better Performance Than Your Current LC Column!". Below the banner is a navigation bar with links for "Easy Column Selection", "Core-Shell Technology", "Technical Resources", "Application Search", "Ask Ball-Guy", "Core Club", and "Order Now". The main content area is titled "Thrive with Kinetex Core-shell Technology" and includes a paragraph about the benefits of core-shell technology, a "FREE ON-SITE Demonstration or Seminar" button, and a "See What Core-shell Technology Can do for Your Lab" link. Below this is a section titled "Complete Scalable Solutions from UHPLC to HPLC to PREP LC" featuring four product images with their particle sizes: 1.3 μm, 1.7 μm, 2.6 μm, and 5 μm. Each image includes a "NEW" badge and a "Buy Now" button. At the bottom of the page are three promotional boxes: "What's New!" with a "Check It Out" button, "Save Now!" with a "Buy Now" button, and "Product Videos" with a "Watch Now" button.

The
Power
of
Core-Shell
is
in
Your
Hands!

IL69770315_W

Phenomenex products are available worldwide

www.phenomenex.com/mysupport

QUALITY MANAGEMENT SYSTEM
CERTIFIED BY DNV

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