

# A Little Joy for Your Corner of the Lab.

Great savings  
on LC, GC, and  
sample prep  
products!



**phenomenex**<sup>®</sup>  
...breaking with tradition<sup>SM</sup>

[www.phenomenex.com/labbuds](http://www.phenomenex.com/labbuds)

# GROW YOUR MIND.

Expand your technical knowledge with our library of on-demand webinars!



Tune into a myriad of informative webinars and learn about trending industry focused solutions, new products to help boost your analysis, and troubleshooting tips!

[www.phenomenex.com/webinars](http://www.phenomenex.com/webinars)

Plant the seeds of your intention & watch it grow!

On-Demand Webinar



# YOUR SAVINGS HAVE BLOSSOMED!

Save on Phenomenex LC and GC columns, and sample preparation products

## HPLC/UHPLC COLUMN DEALS

**10% OFF** any Gemini, Luna, Luna Omega, Kinetex, or Synergi HPLC/UHPLC columns **OR** **15% OFF** 2 or more Gemini, Luna, Luna Omega, Kinetex, or Synergi HPLC/UHPLC columns



Limit 5 columns per company. Excludes guard columns. Offer applies to columns with  $\leq 4.6$  mm ID, and only valid for part numbers purchased for the first time.

## HPLC/UHPLC COLUMN DEALS

**15% OFF** any Aeris, bioZen, Rezex, or Lux LC columns **OR** **20% OFF** 2 or more Aeris, bioZen, Rezex, or Lux LC columns



Limit 5 columns per company. Excludes guard columns. Offer applies to columns with  $\leq 7.8$  mm ID  $\leq 10$  mm ID for Rezex, and only valid for part numbers purchased for the first time.

## GC COLUMN DEALS

**10% OFF** any Zebtron™ GC columns **OR** **15% OFF** 2 or more Zebtron GC columns



Limit 5 columns per company. Excludes custom columns and guard columns. Offer only valid for part numbers purchased for the first time.

## SAMPLE PREPARATION

Buy 2 Boxes, Get 1 Box **FREE**

(must be of equal or lesser value)



Limit 10 boxes per company. Applies to  $\beta$ -Gone™, Novum™, Phree™, roQ™, Strata®, Strata DE, Strata-X, and Strata-X PRO products. Offer only valid for part numbers purchased for the first time.

## HPLC/UHPLC COLUMN PROTECTION

**20% OFF** Your purchase of 2 or more packs of SecurityGuard™ products



Limit 6 products at discounted rate, per company. Offer only valid for part numbers purchased for the first time.

## ACCESSORIES

Contact Your Local Phenomenex Representative!

Offer applies to Phenex™ syringe filters and Verex™ vials.



Limit 5 per company. Offer applies to Phenex and Verex products. Offer only valid for part numbers purchased for the first time.

Offers expire: **April 30, 2021** | Mention offer code: **CATINF21**

Start chatting today at:  
[www.phenomenex.com/chat](http://www.phenomenex.com/chat)

Terms and Conditions  
Offers valid only for addressee. Offers must be specified at time of order and cannot be combined with any other discount or promotional item. All discounts valid for Phenomenex manufactured products only. Offers exclude accessories, columns > 10 mm ID, custom columns, and guard columns. All offers exclude Hypercarb®, PLgel™, PLRP-S™, Spherisorb®, TSKgel®, Ultron®, and ZORBAX® columns. Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

We Get by With a Little Help From Our Friends

Don't know what you need?  
Have technical questions?  
We're here to help!



# LC PRODUCT SPROUTLIGHT

## bioZen LC Columns for Biologics

Maximize sensitivity and selectivity for all your biologics analyses. With the unparalleled inertness of titanium column hardware to minimize priming, three particle platforms for optimal versatility, and eight particle chemistries, bioZen UHPLC/HPLC columns are seamlessly designed to provide enhanced characterization over a wide range of techniques.

- Peptide Mapping (RP)
- Aggregate Analysis (SEC)
- Charge Variant Analysis (IEX)
- Glycan Analysis (HILIC)
- Peptide Quantitation (RP)
- Intact and Fragment Analysis (RP)
- Intact Mass (RP)
- Drug Antibody Ratio (RP)

[www.phenomenex.com/bioZen](http://www.phenomenex.com/bioZen)



## Lux LC Columns for Chiral Separations

Simplify your chiral screening with our collection of Lux amylose and cellulose chiral columns. Lux chiral columns provide a variety of complementary selectivities that allow you to screen under Reversed Phase, Polar Organic, Normal Phase, and SFC conditions.

- Broad enantioselectivity from multiple interaction mechanisms
- 3  $\mu\text{m}$ , 5  $\mu\text{m}$ , 10  $\mu\text{m}$  packed columns for analytical to preparative LC
- Pressure stable up to 300 bar
- High efficiency and loading capacity

[www.phenomenex.com/Lux](http://www.phenomenex.com/Lux)



## Kinetex Core-Shell LC Columns for Small Molecules

Kinetex is a superficially-porous powerhouse for UHPLC, HPLC, and PREP methods. With a vast collection of selectivities, particle and column sizes, Kinetex is a great choice if you need system-to-system and lab-to-lab transferability.

- Scalable from UHPLC to HPLC to Preparative LC
- High efficiencies on any HPLC/UHPLC system
- Transfer easily across all systems and lab-to-lab
- 12 column chemistries for all your methods

[www.phenomenex.com/Kinetex](http://www.phenomenex.com/Kinetex)



## Luna Omega LC Columns for Small Molecules

Luna Omega takes fully-porous technology to astounding levels of performance. An excellent choice for labs using a combination of UHPLC, HPLC, and/or Preparative LC who require a choice of selectivities, and particle and column sizes.

- Scalable from UHPLC to HPLC to Preparative LC
- Rugged and dependable
- Better peak shape through inert foundation
- Several selectivities for method development

[www.phenomenex.com/LunaOmega](http://www.phenomenex.com/LunaOmega)



Discover more LC columns at [www.phenomenex.com/HPLC](http://www.phenomenex.com/HPLC)

# FEATURED GC PRODUCTS

## Zebtron ZB-5MSPLUS™ GC Columns for Versatility

Get next-gen inertness for all your specialty chemical, forensic, toxicology, pharmaceutical, and food testing applications. Experience improved analyte responses and peak shapes with reduced absorption.


- Clean, highly inert fused silica surface
- Specialized deactivation for versatile 5% phenyl-arylene selectivity with improved sensitivity
- Low bleed (MS certified) for high sensitivity GC-MS and GC-MS/MS work



## Zebtron ZB-624PLUS™ GC Columns for Volatiles

The optimal choice for the analysis of volatile compounds. Enjoy low bleed for GC-MS methods, high temperature stability (300/320 C), and enhanced peak shapes due to superior deactivation.

- Excellent column-to-column reproducibility
- Increased sensitivity for boiling point solvents
- Superior deactivation
- Symmetrical peaks for active analytes



Discover more Zebtron phases at [www.phenomenex.com/GC](http://www.phenomenex.com/GC)

# SAMPLE PREPARATION SOLUTIONS

## Strata-X Solid Phase Extraction Products

Remove unwanted contaminants from your samples with this collection of separation mechanisms and separation modes with reliable polymeric sorbent. Choose from reversed phase/ion-exchange/mixed-mode Strata-X sorbents based on your target analytes and achieve longer column lifetimes and reduced MS maintenance.

- Remove proteins, particulates, and de-salt samples
- Tightly retain target analytes
- Clean large sample volumes easily
- Microelution format for limited sample volumes and ultra-concentrated eluates

[www.phenomenex.com/StrataX](http://www.phenomenex.com/StrataX)



## Strata-X PRO Solid Phase Extraction Products

Looking for a faster, cleaner way to extract your samples? Move over traditional SPE, make way for Strata-X PRO! The new revolution in SPE is here and not only needs no method development, but is a great solution for extremely polar analytes.

- Reduce protocol time by at least 40% with 3-step and 2-step SPE
- New Polymeric Sorbent with Matrix Removal Technology
- High recoveries without conditioning or equilibration

[www.phenomenex.com/StrataXPRO](http://www.phenomenex.com/StrataXPRO)



Discover more sample preparation solutions at [www.phenomenex.com/sampleprep](http://www.phenomenex.com/sampleprep)



# A Little Joy for Your Corner of the Lab.



## Australia

t: +61 (0)2-9428-6444  
auiinfo@phenomenex.com

## Austria

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

## Belgium

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

## Canada

t: +1 (800) 543-3681  
info@phenomenex.com

## China

t: +86 400-606-8099  
cninfo@phenomenex.com

## Czech Republic

t: +420 272 017 077  
cz-info@phenomenex.com

## Denmark

t: +45 4824 8048  
nordicinfo@phenomenex.com

## Finland

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

## France

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

## Germany

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

## Hong Kong

t: +852 6012 8162  
hkinfo@phenomenex.com

## India

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

## Ireland

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

## Italy

t: +39 051 6327511  
italiainfo@phenomenex.com

## Luxembourg

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

## Mexico

t: 01-800-844-5226  
tecnicomx@phenomenex.com

## The Netherlands

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

## New Zealand

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

## Norway

t: +47 810 02 005  
nordicinfo@phenomenex.com

## Poland

t: +48 22 104 21 72  
pl-info@phenomenex.com

## Portugal

t: +351 221 450 488  
ptinfo@phenomenex.com

## Singapore

t: +65 800-852-3944  
sginfo@phenomenex.com

## Slovakia

t: +420 272 017 077  
sk-info@phenomenex.com

## Spain

t: +34 91-413-8613  
espinfo@phenomenex.com

## Sweden

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

## Switzerland

t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com

## Taiwan

t: +886 (0) 0801-49-1246  
twinfo@phenomenex.com

## Thailand

t: +66 (0) 2 566 0287  
thaiinfo@phenomenex.com

## United Kingdom

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

## USA

t: +1 (310) 212-0555  
info@phenomenex.com

## 🌐 All other countries/regions

### Corporate Office USA

t: +1 (310) 212-0555  
info@phenomenex.com



## www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

### Terms and Conditions

Offers valid for intended recipient only. Offers must be specified at time of order and may not be combined with any other discount or promotional item. All offers valid for Phenomenex manufactured products and part numbers purchased for the first time only. Offer only valid for part numbers purchased for the first time. Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

### Trademarks

Gemini, Kinetex, Luna, Lux, and Strata are registered trademarks of Phenomenex. Aeris, bioZen, β-Gone, Novum, Phenogel, Phree, Rezex, roQ, SecurityGuard, SecurityLINK, Synergi, Yarra, Zebron, and 5MSPLUS are trademarks of Phenomenex. SelectScience and Scientists' Choice Awards are registered trademarks of SelectScience Ltd.

Gemini and Kinetex EVO are patented by Phenomenex. U.S. Patent Nos. 7,563,367 and 8,658,038 and foreign counterparts.

Strata-X is patented by Phenomenex. U.S. Patent No. 7,119,145.

SecurityGuard is patented by Phenomenex. U.S. Patent No. 6,162,362.

**CAUTION: this patent only applies to the analytical-sized guard cartridge holder, and does not apply to SemiPrep, PREP or ULTRA holders, or to any cartridges.**

Novum is patent pending.

**FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.**

© 2021 Phenomenex, Inc. All rights reserved.